MR6401 CPR Sensor

Operator's Manual

Shanghai International Holding Corp. GmbH(Europe)
Eiffestraße 80, 20537 Hamburg Germany

mindray

| | |
|------------------------|----|
| English | 1 |
| Polski | 9 |
| Deutsch | 17 |
| Français | 25 |
| Nederlands | 33 |
| Brasil Português | 41 |
| Türkçe | 49 |
| Español..... | 57 |
| Italiano | 65 |
| 中文..... | 73 |

mindray is the registered trademark owned by Mindray in China.


©2016-2018 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

Product Introduction

The CPR sensor is intended to provide real-time CPR feedback for patients at least 8 years old or above 25kg weight. It displays the depth and rate of chest compression as well as interruption time. The CPR sensor should be used in healthcare facilities by clinical professionals trained in CPR and use of the device.

The CPR sensor can be configured with or without a battery. It consist of the following accessories:

- MR6801 CPR sensor cable
- MR6802 CPR sensor cable
- MR6921 CPR adhesive tape

Safety Information

WARNING

- Do not use the CPR sensor on patients under 8 years old or below 25kg weight.
 - The CPR sensor is used for single patient at a time.
 - The CPR sensor can be use together with Mindray defibrillator and Mindray monitor supporting the CPR function.
 - When the CPR sensor is used together with a defibrillator, make sure to follow the defibrillator manufacturer's instructions. Stop compressions, remove hands from the CPR sensor and remain clear of all patient contact during defibrillation or when otherwise required, in accordance with a proper defibrillation protocol.
 - When performing CPR on a patient lying on a mattress, a backboard must be used to limit the amount of compressed depth which is absorbed by the mattress. Depending on characteristics of the mattress, backboard and patient, the compensation depth does not guarantee that the patient chest is compressed by 50 mm.
 - The battery indicator lighting in red indicates a low battery. To ensure the CPR sensor is ready for operation in an emergency treatment, you should charge the battery in time. Do not interrupt CPR when the battery indicator lighting in red during the emergency treatment.
 - Do not interrupt CPR in any case, such as low battery or self-test error. Continue CPR without feedback from the CPR sensor.
 - The CPR sensor can bear at most 80kg applied force. Applying force over this limit may result in sensor damage.
 - The CPR sensor is not intended for use in a moving environment, such as an ambulance. If used during patient transport, the CPR sensor may provide inaccurate feedback. If CPR is indicated in a moving environment, do not rely on the CPR sensor depth feedback during such conditions.
 - Do not use the CPR sensor in conjunction with any mechanical or automated compression device.
 - Do not use the CPR sensor on top of defibrillation pads.
 - Cleaning or disinfection should be performed on the CPR sensor housings after each use.
 - If there are signs of fluid entry, stop using the CPR sensor immediately and contact your service personnel.
 - Do not open the CPR sensor housings. All servicing and future upgrades must be carried out by trained and authorized personnel.
 - Replacing the battery if the CPR sensor is configured with is not allowed. It must be replaced by your service personnel.
 - Properly performed CPR may result in fracturing of the patient's ribs, external chest wall bruising or abrasion.
 - If there are externally visible damages on the CPR sensor or the sensor cable, stop using the CPR sensor immediately and contact your service personnel.
-
-

CAUTION

- **Magnetic and electrical fields are capable of interfering with the proper performance of the CPR sensor. For this reason make sure that all external devices operated in the vicinity of the CPR sensor comply with the relevant EMC requirements. Mobile phone, X-ray equipment or MRI devices are a possible source of interference as they may emit higher levels of electromagnetic radiation.**
- **Always carry the CPR sensor properly to avoid damage caused by drop, impact, strong vibration or other mechanical force.**
- **Dry the CPR sensor immediately in case of rain or water spray.**
- **At the end of its service life, the CPR sensor, as well as CPR sensor cable, must be disposed of in compliance with the guidelines regulating the disposal of such products. If you have any questions concerning disposal, please contact us.**
- **Dispose of the package material as per the applicable waste control regulations. Keep it out of children's reach.**

Indicators

| CPR Sensor Status | Compression Depth/ Recoil Indicator | Compression Rate Indicator/ Inactivity Timer |
|---|--|--|
| Power off | Off | Off |
| Sleep (only for CPR sensor configured with a battery) | Off | Off |
| Normal use | indicates the current compression depth, <ul style="list-style-type: none"> ■ Yellow area: 0 to 5 cm ■ Green area: 5 to 6 cm ■ Red area: 6 to 8 cm | indicates the current compression rate |
| Inactivity | the first indicator lights in blue | displays interruption time in seconds since the last compression |
| Self-test error | the first indicator lights in blue | displays EXX (XX represents code 01 to 99) |

Battery

The battery indicators indicate battery status as follows:

- Green: the battery power is normal.
- Yellow: the battery is being charged.
- Red: the battery power is low.
- Off: the battery is depleted or the CPR sensor is not configured with a battery.

When the CPR sensor configured with a battery is in the sleeping status, the battery indicator blinks once every 5 seconds in normal case.

When the battery is low, you should charge the battery by connecting the CPR sensor to a Mindray defibrillator or monitor supporting the CPR function.

WARNING

- **Check the battery indicator regularly.**
 - **The battery indicator lighting in red indicates a low battery. To ensure the CPR sensor is ready for operation in an emergency treatment, you should charge the battery in time. Do not interrupt CPR when the battery indicator lighting in red during the emergency treatment.**
-
-

To store the CPR sensor with a battery for an extended period of time, place the CPR sensor at a temperature of 25°C±5°C and a relative humidity (non-condensing) between 50% and 60%.


For the CPR sensor configured with a battery, you should perform a complete charge-discharge cycle every three months. Make sure the battery is fully charged before the next use.

We recommend replacing the battery every two years.

Basic Operations

1. Apply a new CPR adhesive tape.

Open the package of CPR adhesive tapes and peel one to expose it on the adhesive mount area of the CPR sensor with its adhesive surface.


2. Wake up the CPR sensor.

Only the CPR sensor configured with a battery can enter the sleeping status. In this case, you can wake up the CPR sensor by shaking and compressing it, or connecting it to a Mindray defibrillator or monitor supporting the CPR function.

3. Attach the CPR sensor.

- a. Prepare the patient skin, shave hair from skin if necessary. If the patient's chest is contaminated with blood, clean it completely.
- b. Peel off the liner from the CPR adhesive tape.
- c. Position the CPR sensor so that the compression area is in the middle of the chest, on the lower half of the sternum.

WARNING

- **Do not apply the CPR sensor to an open wound or recent incision site.**
-
-

Cleaning/Disinfection

WARNING

- **Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).**
- **When cleaning or disinfecting the CPR sensor, avoid the CPR sensor connector.**

You should clean or disinfect the CPR sensor as per your hospital's regulations after each use. Sterilization is not allowed for the CPR sensor.

Recommended disinfecting agents are:

- Ethanol (70%)
- Isopropyl alcohol (70%)
- Sodium hypochlorite bleach (10%, Sodium hypochlorite)
- Super Sani-cloth® Germicidal Disposable Wipe (Quaternary ammonium chloride 0.5%, Isopropyl alcohol 55%)
- Hydrogen peroxide (3%)
- Rely+On™ Virkon®
- 1-Propanol (50%)

Troubleshooting

If you encounter problems when using the CPR sensor, check the table below before requesting for services. If the problem persists after you have taken corrective actions, contact your service personnel.

| Problem | Corrective Actions |
|--|---|
| The compression depth/recoil indicator or compression rate indicator/inactivity timer is off. | If your CPR sensor is configured with a battery, check that it is not in the sleeping status. |
| The battery indicator is off. | If it is the CPR sensor configured with a battery, check that its battery power is normal. If the battery is depleted, charge the battery. |
| The compression rate indicator/inactivity timer displays EXX (XX represents code 01 to 99). | Contact your service personnel. |
| The CPR sensor does not adhere to patient's chest. | <ul style="list-style-type: none"> ■ Check that the CPR adhesive tape is not expired. ■ Check that the patient's chest is clean and dry. <p>If the problem persists, administrate CPR using the CPR sensor even if it does not adhere to the patient's chest.</p> |
| You are unfamiliar with the CPR sensor, or doubts there is a problem with the CPR sensor. | <ul style="list-style-type: none"> ■ Ignore feedback from the CPR sensor and continue with CPR. ■ At an appropriate time after the CPR event, check that the CPR adhesive tape is not soiling or damaged. Apply a new CPR adhesive tape and let a trained operator perform chest compressions with the CPR sensor on a CPR manikin. If the CPR sensor does not function as expected, remove it from service and contact your service personnel. |

| Problem | Corrective Actions |
|--|---|
| The sensor cable is not securely inserted, or it easily falls off. | Plug the sensor cable into the connector and fasten it with the cable retainer. |

Product Specifications

WARNING

- **The CPR sensor may not meet the performance specifications if stored or used outside the specified temperature and humidity ranges.**

| | | | |
|---|---|---|----------------------------|
| Safety Classifications (according to IEC 60601-1) | Degree of protection against electrical shock | Type CF defibrillation proof | |
| | Type of protection against electrical shock | <ul style="list-style-type: none"> ■ Used alone: Class II, equipment energized from an internal electrical power source. ■ Used together with the defibrillator or monitor: Class II, equipment energized from a specific or internal electrical power source. | |
| | Degree of protection against harmful ingress of water | IP55 | |
| | Drop height | 1.5 m | |
| | Degree of protection against hazards of explosion | The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air with oxygen or nitrous oxide. | |
| | Mode of operation | Continuous | |
| Environmental Specifications | | Operating conditions | Storage conditions |
| | Temperature | 0°C to 50°C | -30°C to 70°C |
| | Relative humidity | 10% to 95%, non-condensing | 10% to 95%, non-condensing |
| | Barometric | 57.0kPa to 106.2kPa | 57.0kPa to 106.2kPa |
| Power Specifications | Input voltage | 5 to 12V  | |
| | Power consumption | Not applicable | |
| Battery Specifications (optional) | Capacity | ≥230 mAh | |
| | Voltage | 3.8 V | |
| | Run time | <ul style="list-style-type: none"> ■ At least 3 hours when the CPR sensor works continuously at a temperature of 25°C±5°C with a fully-charged battery. ■ At least 30 days when the CPR sensor is in the sleeping status. ■ At least 1 hour after the battery indicator is in red light. | |
| | Charge time | At most 1.5 hours to a full charge when a depleted battery is charged at a temperature of 25°C±5°C. | |

| | | | |
|----------------------------|---------------------|--|---|
| Physical Specifications | Sensor size | < 150mm × 63mm × 20m | |
| | Compression area | < 92mm × 53mm × 19mm | |
| | Adhesive mount area | ≤ 98mm × 45mm | |
| | Sensor weight | ≤ 220g (including battery, without CPR sensor cable) | |
| Measurement Specifications | Compression depth | Measurement range | 0.0 to 8.0 cm |
| | | Effective range | 1.5 to 8.0 cm |
| | | Accuracy | ±0.5 cm or ±10%, whichever is greater |
| | | Resolution | 0.1 cm |
| | | Refreshing rate | ≥ 0.5 Hz |
| | Compression rate | Measurement range | 40 to 160 cpm (compressions per minute) |
| | | Effective range | 40 to 160 cpm (compressions per minute) |
| | | Accuracy | ±2 cpm (compressions per minute) |
| | | Resolution | 1 cpm |
| | | Refreshing rate | ≥ 0.5 Hz |
| | Interruption time | Measurement range | 0 to 300 s |
| | | Effective range | 0 to 300 s |
| | | Resolution | 1 s |
| Refreshing rate | | ≥ 0.5 Hz | |

EMC

The device meets the requirements of IEC 60601-1-2. For more information, see the operator's manual of the defibrillator or monitor used together with this device.

Symbol

| | | | |
|---|---|---|---|
|  | Refer to instruction manual/booklet |  | Date of manufacture |
|  | Serial number | IP55 | Dust-protected Protected against water jets |
|  | DEFIBRILLATION-PROOF TYPE CF APPLIED PART | EC REP | Authorized representative in the European Community |
|  | Not made with natural rubber latex. | LOT | Batch code |
|  | Manufacturer |  | Attention, consult ACCOMPANYING DOCUMENTS |

MR6401 CPR Sensor Operator's Manual

| | |
|---|--|
|  | Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive). |
|  | Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive). |
|  | Comply with the requirements of Directive 2012/19/EU Waste Electrical & Electronic Equipment |

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