# Multifunction Electrode Pads

## Instructions for Use

|                 |    |
|-----------------|----|
| English .....   | 1  |
| Deutsch .....   | 7  |
| Français .....  | 13 |
| Nederlands..... | 19 |
| Italiano .....  | 25 |
| Español .....   | 31 |
| Português ..... | 37 |
| Türkçe .....    | 43 |
| Polski .....    | 49 |
| Русский .....   | 55 |
| Česky .....     | 61 |
| Magyar.....     | 67 |
| 中文 .....        | 73 |

**mindray**

## Statement

**SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.** (hereinafter called Mindray) owns the intellectual property rights to this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- the product is used in accordance with the instructions for use.
- this product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.

In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief dialescription of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

**mindray.**  **MINDRAY** and  are the registered trademarks or trademarks owned by Mindray in China.

© 2009-2019 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

## Product Introduction

### Intended Use

MR61/MR63 multifunction electrode pads are pre-gelled, self adhesive electrodes intended for ECG monitoring, noninvasive pacing, manual defibrillation and synchronized cardioversion for pediatric patients with weight less than 25kg.

This product must only be used by trained clinical professionals and with Mindray-specified defibrillator/monitor. The operator shall be liable for understanding this instructions for use or consulting with us to confirm the compatibility of this product and the defibrillator/monitor. Otherwise, product performance may be compromised or personal injury may occur. This product can be used for up to 25 defibrillation shocks.

### Warnings

- Defibrillation energy shall not be greater than 100J.
- Do not crush, bend or fold electrode pads, or place them under heavy objects.
- Check expiration date on package. Do not use if beyond the expiration date.
- Check pads package for integrity before use. Remove from use if there is sign of damage.
- Do not open pads package until immediately prior to use.
- Check the electrodes before use. Remove them from use in case of damage, deterioration or visible contamination.
- Do not use if gel is dry or separated from the metal backing.
- Do not use additional conductive gel and paste on electrodes.
- Do not apply electrode pads over broken skin, if possible.
- Keep electrode pads clear of any other electrodes or metal parts in contact with the patient during defibrillation. Change electrode pads after 25 defibrillation shocks.

## Warnings

- Check patient skin at application site after 30 minutes of continuous pacing. Replace electrode pads after continuous pacing for 1 hours.
- Electrode pads are not repositionable. Replace with new pads if repositioning is required.
- Change electrode pads every 8 hours.
- Do not clean electrodes with alcohol or other cleaning agents.
- Do not immerse cables or electrode pads.
- The product is non-sterile. Do not sterilize.
- For single patient use only. Discard after use.


## Applying the Electrode Pads

Be sure to understand above precautions and warnings before using the multifunction electrode pads.

1. Select placement sites on body of patient according to defibrillator/monitor manufacturer's operating manual. Figure A and Figure B shows two standard opations.
2. Clean and dry the skin of application site. Remove excessive hair if necessary.
3. Remove pads from package. Separate leadwires. Peel the protective cover of electrode.
4. Apply electrodes as directed by Figure A or Figure B. Smooth the electrode from the center outward to edges with fingertips to ensure that there are no air pockets between the gel and patient's skin.
5. Connect the electrodes to therapy cable and then connect the cable to defibrillator/monitor.
6. Follow the defibrillator/monitor manufacturer's operating manual to perform patient therapy or ECG monitoring. Replace electrode pads after 25 defibrillation shocks, 1 hour of continuous pacing or 8 hours on skin.


MULTIFUNCTION ELECTRODE PADS INSTRUCTIONS FOR USE

A







Anterior-lateral placement

B

Anterior-posterior placement


## SPEC/Symbol

### Storage

Storage Temperature: 0°C to 50°C (32°F to 122°F).

### Symbols

|   |   |  |   |
|---|---|--|---|
|  | Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive). |                   | Refer to instruction manual/ booklet                  |
|  | Do not re-use   |  <b>1d</b><br>MAX | Use within 24 hours if unpacked                       |
|  | Non-sterile   |                   | Keep dry<br>Keep away from sunlight                   |
|  | Manufacturer  |                   | Temperature limit<br>Do not use if package is damaged |
|  | Date of manufacture   |                   | Use-by date   |
|  | Batch code  |                   | Caution   |

### Rx Only

Federal law (USA) restricts this device to sale by or on the order of a physician.

**This page is intentionally left blank.**

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