Multifunction Electrode Pads

Instructions for Use

English	1
Deutsch	7
Français	13
Nederlands.....	19
Italiano	25
Español	31
Português	37
Türkçe	43
Polski	49
Русский	55
Česky	61
Magyar	67
中文	73

mindray

Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety.

Mindray is responsible for the effects on safety, reliability and performance of this product, only if:

- the product is used in accordance with the instructions for use.
- this product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.

In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief dialescription of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

mindray.  **MINDRAY** and  are the registered trademarks or trademarks owned by Mindray in China.

© 2009-2019 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

Product Introduction

Intended Use

MR60/MR62 multifunction electrode pads are pre-gelled, self adhesive electrodes intended for ECG monitoring, noninvasive pacing, automated external defibrillation, manual defibrillation and synchronized cardioversion for adult patients.

This product must only be used by trained clinical professionals and with Mindray-specified defibrillator/monitor. The operator shall be liable for understanding this instructions for use or consulting with us to confirm the compatibility of this product and the defibrillator/monitor. Otherwise, product performance may be compromised or personal injury may occur.

Warnings

- Do not crush, bend or fold the pads, or place them under heavy objects.
- Check expiration date on package. Do not use if beyond the expiration date.
- Check pads package for integrity before use. Remove from use if there is sign of damage.
- Do not open pads package until immediately prior to use.
- Check the pads before use. Remove them from use in case of damage, deterioration or visible contamination.
- Do not use if gel is dry or separated from the metal backing.
- Do not use additional conductive gel and paste on pads.
- Do not apply pads over broken skin, if possible.
- Keep pads clear of any other electrodes or metal parts in contact with the patient during defibrillation.
- Check patient skin at application site every 30 minutes of continuous pacing. Replace pads after continuous pacing for 8 hours.

Warnings

- Electrode pads are not repositionable. Replace with new pads if repositioning is required.
- Change pads after 50 defibrillation shocks.
- Change pads every 24 hours.
- Do not clean pads with alcohol or other cleaning agents.
- Do not immerse cables or electrode pads.
- The product is non-sterile. Do not sterilize.
- For single patient use only. Discard after use.

Applying the Electrode Pads


Be sure to understand above precautions and warnings before using the multifunction electrode pads.

1. Select placement sites on body of patient according to defibrillator/monitor manufacturer's operating manual. Figure A and Figure B shows two standard opations.
2. Clean and dry the skin of application site. Remove excessive hair if necessary.
3. Remove electrode pads from package. Separate leadwires. Peel the protective cover of electrode.
4. Apply electrode pads as directed by Figure A or Figure B. Smooth the electrode from the center outward to edges with fingertips to ensure that there are no air pockets between the gel and patient's skin.
5. Connect the electrodes to therapy cable and then connect the cable to defibrillator/monitor.
6. Follow the defibrillator/monitor manufacturer's operating manual to perform patient therapy or ECG monitoring.. Replace electrode pads after 50 defibrillation shocks, 8 hours of continuou pacing or 24 hours on skin.







MULTIFUNCTION ELECTRODE PADS INSTRUCTIONS FOR USE

Anterior-lateral placement

Anterior-posterior placement

SPEC/Symbol

Storage

Storage Temperature: 0°C to 50°C (32°F to 122°F).

Symbols

Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive).

Refer to instruction manual/ booklet

Do not re-use

Use within 24 hours if unpacked

Non-sterile

Keep dry

Keep away from sunlight

Manufacturer

Temperature limit
Do not use if package is damaged

Date of manufacture

Use-by date

Batch code

Caution

Rx Only

Federal law (USA) restricts this device to sale by or on the order of a physician.

This page is intentionally left blank.

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