

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200**: BS-200, BS-220;
2. **BS-200E**: BS-200E, BS-220E;
3. **BS-230**: BS-230, BS-240;
4. **BS-240E**: BS240E, BS240Pro;
5. **BS-300**: BS-300, BS-320;
6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-");
7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E (V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E (V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

8. **BS-380**: BS-380, BS-390;
9. **BS-400**: BS-400, BS-420;
10. **BS-430**: BS-430, BS-450, BS-460;
11. **BS-480**: BS-480, BS-490;
12. **BS-600**: BS-600, BS-620;
13. **BS-600M**: BS-600M, BS-620M;
14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;
16. **BS-2800M**: BS-2600M, BS-2800M.

17. **S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 145222002 : 2023-10-31

Abbreviated name			CHE		
Model	Unit	Calibration Value ¹⁷	Model	Unit	Calibration Value ¹⁷
BS-200 ¹	U/L	5200	BS-400 ⁹	U/L	5200
	µkat/L	86.8		µkat/L	86.8
BS-200E ²	U/L	5200	BS-430 ¹⁰	U/L	5200
	µkat/L	86.8		µkat/L	86.8
BS-230 ³	U/L	5200	BS-480 ¹¹	U/L	5200
	µkat/L	86.8		µkat/L	86.8
BS-240E ⁴	U/L	5200	BS-600 ¹²	U/L	5200
	µkat/L	86.8		µkat/L	86.8
BS-300 ⁵	U/L	5200	BS-600M ¹³	U/L	5243
	µkat/L	86.8		µkat/L	87.6
BS-330E ⁶	U/L	5200	BS-800 ¹⁴	U/L	5200
	µkat/L	86.8		µkat/L	86.8
BS-360E ⁷	U/L	5200	BS-2000 ¹⁵	U/L	5247
	µkat/L	86.8		µkat/L	87.6
BS-380 ⁸	U/L	5200	BS-2800M ¹⁶	U/L	5243
	µkat/L	86.8		µkat/L	87.6

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-200**: BS-200, BS-220;

2. **BS-200E**: BS-200E, BS-220E;

3. **BS-230**: BS-230, BS-240;

4. **BS-240E**: BS240E, BS240Pro;

5. **BS-300**: BS-300, BS-320;

6. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),
BS-350E (Serial Number starts with "XS-");

7. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E (V35.00) (Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E (V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

8. **BS-380**: BS-380, BS-390;

9. **BS-400**: BS-400, BS-420;

10. **BS-430**: BS-430, BS-450, BS-460;

11. **BS-480**: BS-480, BS-490;

12. **BS-600**: BS-600, BS-620;

13. **BS-600M**: BS-600M, BS-620M;

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 145222002 : 2023-10-31

CHE													
Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value ± 3SD)		1 SD		
BS-200	1 U/L	5232	4180	—	6284	351	BS-400	9 U/L	5256	4200	—	6312	352
	µkat/L	87.4	69.8	—	104.9	5.9		µkat/L	87.8	70.1	—	105.4	5.9
BS-200E	2 U/L	5122	4092	—	6152	343	BS-430	10 U/L	5266	4208	—	6324	353
	µkat/L	85.5	68.3	—	102.7	5.7		µkat/L	87.9	70.3	—	105.6	5.9
BS-230	3 U/L	5217	4168	—	6266	350	BS-480	11 U/L	5297	4232	—	6362	355
	µkat/L	87.1	69.6	—	104.6	5.8		µkat/L	88.5	70.7	—	106.2	5.9
BS-240E	4 U/L	5226	4176	—	6276	350	BS-600	12 U/L	5266	4208	—	6324	353
	µkat/L	87.3	69.7	—	104.8	5.8		µkat/L	87.9	70.3	—	105.6	5.9
BS-300	5 U/L	5173	4133	—	6213	347	BS-600M	13 U/L	5255	4199	—	6311	352
	µkat/L	86.4	69.0	—	103.8	5.8		µkat/L	87.8	70.1	—	105.4	5.9
BS-330E	6 U/L	5122	4092	—	6152	343	BS-800	14 U/L	5266	4208	—	6324	353
	µkat/L	85.5	68.3	—	102.7	5.7		µkat/L	87.9	70.3	—	105.6	5.9
BS-360E	7 U/L	5279	4218	—	6340	354	BS-2000	15 U/L	5271	4212	—	6330	353
	µkat/L	88.2	70.4	—	105.9	5.9		µkat/L	88.0	70.3	—	105.7	5.9
BS-380	8 U/L	5256	4200	—	6312	352	BS-2800M	16 U/L	5255	4199	—	6311	352
	µkat/L	87.8	70.1	—	105.4	5.9		µkat/L	87.8	70.1	—	105.4	5.9

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value ± 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ± 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ± 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ± 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri ± 3SD)

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	colinesterasi
Türkçe	Kolinesteraz