

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

2. **BS-180**: BS-180, BS-190;

11. **BS-380**: BS-380, BS-390;

3. **BS-200**: BS-200, BS-220;

12. **BS-400**: BS-400, BS-420;

4. **BS-200E**: BS-200E, BS-220E;

13. **BS-430**: BS-430, BS-450, BS-460;

5. **BS-230**: BS-230, BS-240;

14. **BS-480**: BS-480, BS-490;

6. **BS-240E**: BS240E, BS240Pro;

15. **BS-600**: BS-600, BS-620;

7. **BS-300**: BS-300, BS-320;

16. **BS-600M**: BS-600M, BS-620M;

8. **BS-330**: BS-330, BS-350;

17. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

9. **BS-330E**: BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-");

18. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

19. **BS-2800M**: BS-2600M, BS-2800M.

20. **S0: 0,9% NaCl, Conc. Of S0=0;**

Русский: S0: 0, 9% NaCl, конц. S0=0;

Português: S0: 0,9% NaCl, Conc. de S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

**LOT** : 145322001

**EXP** : 2023-05-13

Abbreviated name			Fe		
Model	Unit	Calibration Value <sup>20</sup>	Model	Unit	Calibration Value <sup>20</sup>
<b>BS-120</b> <sup>1</sup>	µmol/L	/	<b>BS-380</b> <sup>11</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-180</b> <sup>2</sup>	µmol/L	/	<b>BS-400</b> <sup>12</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-200</b> <sup>3</sup>	µmol/L	/	<b>BS-430</b> <sup>13</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-200E</b> <sup>4</sup>	µmol/L	/	<b>BS-480</b> <sup>14</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-230</b> <sup>5</sup>	µmol/L	/	<b>BS-600</b> <sup>15</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-240E</b> <sup>6</sup>	µmol/L	/	<b>BS-600M</b> <sup>16</sup>	µmol/L	41.2
	mg/L	/		mg/L	2.31
<b>BS-300</b> <sup>7</sup>	µmol/L	/	<b>BS-800</b> <sup>17</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-330</b> <sup>8</sup>	µmol/L	/	<b>BS-2000</b> <sup>18</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-330E</b> <sup>9</sup>	µmol/L	/	<b>BS-2800M</b> <sup>19</sup>	µmol/L	/
	mg/L	/		mg/L	/
<b>BS-360E</b> <sup>10</sup>	µmol/L	/			
	mg/L	/			

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

		Fe
English		Iron
Русский		Железо
Português		Ferro
Español		Hierro
Italiano		Ferro
Türkçe		Demir