

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах. Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma. Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.**BS-230**: BS-230, BS-240;

2.**BS-240E**: BS-240E, BS-240Pro;

3.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

**BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"),**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

4.**BS-380**: BS-380, BS-390;

13. **Cal0: 0.9% NaCl, Conc. Of Cal0=0.**

Português: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Italiano: Cal0: 0,9% NaCl, conc. di Cal0=0.

5.**BS-400**: BS-400, BS-420;

6.**BS-430**: BS-430, BS-450, BS-460;

7.**BS-480**: BS-480, BS-490;

8.**BS-600**: BS-600, BS-620;

9.**BS-600M**: BS-600M, BS-620M;

10.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

11.**BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

12.**BS-2800M**: BS-2600M, BS-2800M.

Русский: Cal0: 0,9% NaCl, конц. Cal0=0.

Español: Cal0: 0,9% NaCl, Conc. de Cal0=0.

Türkçe: Cal0: %0,9 NaCl, Cal0 Kons.=0.

**LOT** : 146723005

**EXP** : 2024-08-04

| Abbreviated name              |             | UIBC                            |       |
|-------------------------------|-------------|---------------------------------|-------|
| Model                         | Level       | Calibration Value <sup>13</sup> |       |
|                               |             | µmol/L                          | µg/dL |
| <b>BS-230</b> <sup>1</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-240E</b> <sup>2</sup>   | <b>Cal1</b> | 117                             | 654   |
| <b>BS-360E</b> <sup>3</sup>   | <b>Cal1</b> | 117                             | 654   |
| <b>BS-380</b> <sup>4</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-400</b> <sup>5</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-430</b> <sup>6</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-480</b> <sup>7</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-600</b> <sup>8</sup>    | <b>Cal1</b> | 117                             | 654   |
| <b>BS-600M</b> <sup>9</sup>   | <b>Cal1</b> | 112                             | 626   |
| <b>BS-800</b> <sup>10</sup>   | <b>Cal1</b> | 117                             | 654   |
| <b>BS-2000</b> <sup>11</sup>  | <b>Cal1</b> | 114                             | 637   |
| <b>BS-2800M</b> <sup>12</sup> | <b>Cal1</b> | 112                             | 626   |

| English   | Abbreviated name         | Calibration Value      | Model   | Level   |
|-----------|--------------------------|------------------------|---------|---------|
| Русский   | сокращенное наименование | Эталонное значение     | модель  | Уровень |
| Português | Nome abreviado           | Valor de calibração    | Modelo  | Nível   |
| Español   | nombre abreviado         | Valor de calibración   | modelo  | Nivel   |
| Italiano  | abbreviazione            | Valore di calibrazione | modelli | Livello |
| Türkçe    | kısaltılmış ad           | Kalibrasyon Değeri     | model   | Düzey   |

| UIBC      |  |
|-----------|--|
| English   | Unsaturated Iron Binding Capacity          |
| Русский   | ненасыщенная железосвязывающая способность |
| Português | Capacidade de ligação de ferro insaturado  |
| Español   | Capacidad de unión de hierro no saturado   |
| Italiano  | Capacità di legame del ferro insaturo      |
| Türkçe    | Unsature Demir Bağlama Kapasitesi          |