

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a r

Italiano: la dati di ogni gruppo è la ste

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-330**: BS-330, BS-350;

7. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,
BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00"),
BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

15. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

16. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 147922004

: 2023-09-06

| Abbreviated name | Model | Unit | Calibration Value | Model | Unit | Calibration Value |
|------------------|-----------------------------|------|-------------------|-------------------------------|------|-------------------|
| TPUC | BS-120 ¹ | mg/L | 996 | BS-380 ⁹ | mg/L | 1042 |
| | BS-200 ² | mg/L | 1011 | BS-400 ¹⁰ | mg/L | 1042 |
| | BS-200E ³ | mg/L | 1042 | BS-430 ¹¹ | mg/L | 1042 |
| | BS-230 ⁴ | mg/L | 1000 | BS-480 ¹² | mg/L | 1042 |
| | BS-240E ⁵ | mg/L | 1042 | BS-600 ¹³ | mg/L | 1042 |
| | BS-330 ⁶ | mg/L | 1011 | BS-800 ¹⁴ | mg/L | 1042 |
| | BS-330E ⁷ | mg/L | 1042 | BS-2000 ¹⁵ | mg/L | 1059 |
| | BS-360E ⁸ | mg/L | 1042 | BS-2800M ¹⁶ | mg/L | 1044 |

| English | Abbreviated name | Calibration Value | Unit |
|-----------|--------------------------|------------------------|---------|
| Русский | сокращенное наименование | Эталонное значение | Прибор |
| Português | Nome abreviado | Valor de calibração | Unidade |
| Español | nombre abreviado | Valor de calibración | Unidad |
| Italiano | abbreviazione | Valore di calibrazione | Unità |
| Türkçe | kısaltılmış ad | Kalibrasyon Değeri | Ünite |

| TPUC | |
|-----------|--------------------------------------|
| English | Total Protein in Urine/CSF (TPUC) |
| Русский | Общий белок в моче/СМЖ (TPUC) |
| Português | Proteína total em urina e LCR (PTUL) |
| Español | proteína total en orina/LCR (TPUC) |
| Italiano | Proteine Totali Urinarie/CSF (TPUC) |
| Türkçe | BOS/İdrar Protein (TPUC) |