

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

11.BS-430: BS-430, BS-450, BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

12.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

13.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

14.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-300: BS-300, BS-320;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00)

17.BS-2800M:BS-2600M, BS-2800M.

(Serial Number starts with "W9-" and software version starts with "35.00");

18.0.9% NaCl is needed, Conc. Of 0.9% NaCl =0;

9.BS-380: BS-380, BS-390;

19.Applicable models of the chemistry shall be subject to the

10.BS-400: BS-400, BS-420;

parameter list and instructions.



LOT : 151423010



: 2025-06-06

English	Abbreviated name	Unit	Calibration Rule	Model
Русский	сокращенное наименование	Прибор	Принцип калибровки	модель
Português	Nome abreviado	Unidade	Regra de calibração	Modelo
Español	nombre abreviado	Unidad	Regla de calibración	modelo
Italiano	abbreviazione	Unità	Regola di calibrazione	modelli
Türkçe	kısaltılmış ad	Ünite	Kalibrasyon Kuralı	model

English	Expanded Uncertainty
Русский	Расширенная неопределенность
Português	Incerteza Expandida
Español	Incertidumbre expandida
Italiano	Incertezza estesa
Türkçe	Daha Uzun Süreli Belirsizlik

CK-MB	
English	Creatine Kinase-MB
Русский	MB фракцию креатинкиназы
Português	creatina quinase-MB
Español	creatina quinasa-MB
Italiano	creatina chinasi-MB
Türkçe	Kreatin Kinaz-MB

CK-MB Calibrator



Abbreviated name	Unit	Calibration Rule				Two-point linear	
		Model	Calibrator Value	Expanded Uncertainty	Model	Calibrator Value	Expanded Uncertainty
CK-MB	U/L	BS-120 ¹	119	5	BS-400 ¹⁰	123	5
		BS-200 ²	119	5	BS-430 ¹¹	123	5
		BS-200E ³	123	5	BS-480 ¹²	123	5
		BS-240 ⁴	119	5	BS-600 ¹³	123	5
		BS-240E ⁵	123	5	BS-600M ¹⁴	125	5
		BS-300 ⁶	123	5	BS-800 ¹⁵	123	5
		BS-330E ⁷	123	5	BS-2000 ¹⁶	123	5
		BS-360E ⁸	123	5	BS-2800M ¹⁷	125	5
		BS-380 ⁹	123	5			
	μkat/L	BS-120 ¹	1.99	0.08	BS-400 ¹⁰	2.05	0.08
		BS-200 ²	1.99	0.08	BS-430 ¹¹	2.05	0.08
		BS-200E ³	2.05	0.08	BS-480 ¹²	2.05	0.08
		BS-240 ⁴	1.99	0.08	BS-600 ¹³	2.05	0.08
		BS-240E ⁵	2.05	0.08	BS-600M ¹⁴	2.09	0.08
		BS-300 ⁶	2.05	0.08	BS-800 ¹⁵	2.05	0.08
		BS-330E ⁷	2.05	0.08	BS-2000 ¹⁶	2.05	0.08
		BS-360E ⁸	2.05	0.08	BS-2800M ¹⁷	2.09	0.08
BS-380 ⁹	2.05	0.08					

mindray

CK-MB Calibrator

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 151423010

2025-06-06

