

## Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-120**: BS-120, BS-130;

2.**BS-180**: BS-180, BS-190;

3.**BS-200**: BS-200, BS-220;

4.**BS-200E**: BS-200E, BS-220E;

5.**BS-230**: BS-230, BS-240;

6.**BS-240E**: BS240E, BS240Pro;

7.**BS-300**: BS-300, BS-320;

8.**BS-330**: BS-330, BS-350;

9.**BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

**BS-350E**(Serial Number starts with "XS-");

10.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

**BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"),

**BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

11.**BS-380**: BS-380, BS-390;

12.**BS-400**: BS-400, BS-420;

13.**BS-430**: BS-430, BS-450, BS-460;

14.**BS-480**: BS-480, BS-490;

15.**BS-600**: BS-600, BS-620;

16.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

17.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

18.**BS-2800M**:BS-2600M, BS-2800M.

**LOT** : 155822001

**EXP** : 2023-08-15

Abbreviated name	Model	Unit	Assay Value	Model	Unit	Assay Value
ADA	<b>BS-120</b> <sup>1</sup>	U/L	47.4	<b>BS-360E</b> <sup>10</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-180</b> <sup>2</sup>	U/L	47.4	<b>BS-380</b> <sup>11</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-200</b> <sup>3</sup>	U/L	47.4	<b>BS-400</b> <sup>12</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-200E</b> <sup>4</sup>	U/L	47.4	<b>BS-430</b> <sup>13</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-230</b> <sup>5</sup>	U/L	47.4	<b>BS-480</b> <sup>14</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-240E</b> <sup>6</sup>	U/L	47.4	<b>BS-600</b> <sup>15</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-300</b> <sup>7</sup>	U/L	47.4	<b>BS-800</b> <sup>16</sup>	U/L	47.4
		µkat/L	0.792		µkat/L	0.792
	<b>BS-330</b> <sup>8</sup>	U/L	47.4	<b>BS-2000</b> <sup>17</sup>	U/L	48.2
		µkat/L	0.792		µkat/L	0.805
	<b>BS-330E</b> <sup>9</sup>	U/L	47.4	<b>BS-2800M</b> <sup>18</sup>	U/L	48.4
		µkat/L	0.792		µkat/L	0.808

English	Abbreviated name	Model	Unit	Assay Value
Русский	сокращенное наименование	модель	Прибор	Результат анализа
Português	Nome abreviado	Modelo	Unidade	Valores da análise
Español	nombre abreviado	modelo	Unidad	Valor de ensayo
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio
Türkçe	kisaltılmış ad	model	Ünite	Tayin Değeri

ADA				
English	Adenosine Deaminase			
Русский	является аденозиндезаминаза			
Português	Adenosina deaminase			
Español	adenosina desaminasa			
Italiano	adenosina deamminasi			
Türkçe	Adenozin Deaminazdır			