

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS240E, BS240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-"

and software version starts with "35.00");

14. **S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

5. **BS-380**: BS-380, BS-390;

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 158623002

EXP : 2024-09-29

| Abbreviated name | | CysC II | Calibration Rule | | Spline |
|----------------------|-------|---|-----------------------|-------|---|
| Model | Level | Calibration Value ¹⁴ mg/L | Model | Level | Calibration Value ¹⁴ mg/L |
| BS-200E ¹ | a | 0.590 | BS-480 ⁸ | a | 0.590 |
| | b | 1.16 | | b | 1.16 |
| | c | 2.27 | | c | 2.27 |
| | d | 4.64 | | d | 4.64 |
| | e | 8.62 | | e | 8.62 |
| BS-230 ² | a | 0.590 | BS-600 ⁹ | a | 0.590 |
| | b | 1.16 | | b | 1.16 |
| | c | 2.27 | | c | 2.27 |
| | d | 4.64 | | d | 4.64 |
| | e | 8.62 | | e | 8.62 |
| BS-240E ³ | a | 0.590 | BS-600M ¹⁰ | a | 0.600 |
| | b | 1.16 | | b | 1.17 |
| | c | 2.27 | | c | 2.29 |
| | d | 4.64 | | d | 4.70 |
| | e | 8.62 | | e | 8.82 |
| BS-360E ⁴ | a | 0.590 | BS-800 ¹¹ | a | 0.590 |
| | b | 1.16 | | b | 1.16 |
| | c | 2.27 | | c | 2.27 |
| | d | 4.64 | | d | 4.64 |
| | e | 8.62 | | e | 8.62 |

Cystatin C Calibrator

| Abbreviated name | | CysC II | Calibration Rule | | Spline | | |
|---------------------|-------|---------------------------------|------------------|------------------------|--------|---------------------------------|--|
| Model | Level | Calibration Value ¹⁴ | | Model | Level | Calibration Value ¹⁴ | |
| | | mg/L | | | | mg/L | |
| BS-380 ⁵ | a | 0.590 | | BS-2000 ¹² | a | 0.600 | |
| | b | 1.16 | | | b | 1.20 | |
| | c | 2.27 | | | c | 2.34 | |
| | d | 4.64 | | | d | 4.75 | |
| | e | 8.62 | | | e | 8.78 | |
| BS-400 ⁶ | a | 0.590 | | BS-2800M ¹³ | a | 0.600 | |
| | b | 1.16 | | | b | 1.17 | |
| | c | 2.27 | | | c | 2.29 | |
| | d | 4.64 | | | d | 4.70 | |
| | e | 8.62 | | | e | 8.82 | |
| BS-430 ⁷ | a | 0.590 | | | | | |
| | b | 1.16 | | | | | |
| | c | 2.27 | | | | | |
| | d | 4.64 | | | | | |
| | e | 8.62 | | | | | |

| English | Abbreviated name | Calibration Rule | Model | Level | Calibration Value |
|-----------|--------------------------|------------------------|---------|---------|------------------------|
| Русский | сокращенное наименование | Принцип калибровки | модель | Уровень | Эталонное значение |
| Português | Nome abreviado | Regra de calibração | Modelo | Nível | Valor de calibração |
| Español | nombre abreviado | Regla de calibración | modelo | Nivel | Valor de calibración |
| Italiano | abbreviazione | Regola di calibrazione | modelli | Livello | Valore di calibrazione |
| Türkçe | kısaltılmış ad | Kalibrasyon Kuralı | model | Düzey | Kalibrasyon Değeri |

| CysC | |
|-----------|-------------|
| English | Cystatin C |
| Русский | цистатин C |
| Português | Cistatina C |
| Español | cistatina C |
| Italiano | cistatina C |
| Türkçe | Cystatin C |