

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

1.**BS-120**: BS-120, BS-130;

2.**BS-200**: BS-200, BS-220;

3.**BS-200E**: BS-200E, BS-220E;

4.**BS-230**: BS-230, BS-240;

5.**BS-240E**: BS-240E, BS-240Pro;

6.**BS-300**: BS-300, BS-320;

7.**BS-330E**:BS-330E(Serial Number starts with "XQ-"),

BS-350E(Serial Number starts with "XS-");

18.**S0: 0.9% NaCl, Conc. Of S0=0;**

Português: S0: 0,9% NaCl, Conc. de S0=0;

Italiano: S0: 0,9% NaCl, conc. di S0=0;

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

8.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

9.**BS-380**: BS-380, BS-390;

10.**BS-400**: BS-400, BS-420;

11.**BS-430**: BS-430, BS-450, BS-460;

12.**BS-480**: BS-480, BS-490;

13.**BS-600**: BS-600, BS-620;

14.**BS-600M**: BS-600M, BS-620M;

15.**BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

16.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17.**BS-2800M**: BS-2600M, BS-2800M.

Русский: S0: 0, 9% NaCl, конц. S0=0;

Español: S0: 0,9% NaCl, Conc. de S0=0;

Türkçe: S0: %0,9 NaCl, S0 Kons.=0.

LOT : 044823008

EXP : 2025-01-30

Abbreviated name			LIP		
Model	Unit	Calibration Value ¹⁸	Model	Unit	Calibration Value ¹⁸
BS-120 ¹	U/L	92.5	BS-400 ¹⁰	U/L	95.4
	µkat/L	1.54		µkat/L	1.59
BS-200 ²	U/L	92.5	BS-430 ¹¹	U/L	95.4
	µkat/L	1.54		µkat/L	1.59
BS-200E ³	U/L	95.4	BS-480 ¹²	U/L	95.4
	µkat/L	1.59		µkat/L	1.59
BS-230 ⁴	U/L	92.5	BS-600 ¹³	U/L	95.4
	µkat/L	1.54		µkat/L	1.59
BS-240E ⁵	U/L	95.4	BS-600M ¹⁴	U/L	96.1
	µkat/L	1.59		µkat/L	1.60
BS-300 ⁶	U/L	95.4	BS-800 ¹⁵	U/L	95.4
	µkat/L	1.59		µkat/L	1.59
BS-330E ⁷	U/L	95.4	BS-2000 ¹⁶	U/L	95.6
	µkat/L	1.59		µkat/L	1.60
BS-360E ⁸	U/L	95.4	BS-2800M ¹⁷	U/L	96.1
	µkat/L	1.59		µkat/L	1.60
BS-380 ⁹	U/L	95.4			
	µkat/L	1.59			

English	Abbreviated name	Model	Unit	Calibration Value
Русский	сокращенное наименование	модель	Прибор	Принцип калибровки
Português	Nome abreviado	Modelo	Unidade	Regra de calibração
Español	nombre abreviado	modelo	Unidad	Regla de calibración
Italiano	abbreviazione	modelli	Unità	Regola di calibrazione
Türkçe	kısaltılmış ad	model	Ünite	Kalibrasyon Kuralı

LIP	
English	Lipase
Русский	Липаза
Português	Lipase
Español	Lipasa
Italiano	Lipasi
Türkçe	Lipaz

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-300**: BS-300, BS-320;

7. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

BS-350E (Serial Number starts with "XS-");

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-600M**: BS-600M, BS-620M;

15. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

16. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 044823008

EXP : 2025-01-30

Abbreviated name						LIP							
Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
BS-120 ¹	U/L	47.9	38.3	—	57.5	3.2	BS-400 ¹⁰	U/L	45.4	36.4	—	54.4	3.0
	μ kat/L	0.800	0.640	—	0.960	0.053		μ kat/L	0.758	0.608	—	0.908	0.050
BS-200 ²	U/L	43.2	34.5	—	51.9	2.9	BS-430 ¹¹	U/L	46.1	36.8	—	55.4	3.1
	μ kat/L	0.721	0.576	—	0.867	0.048		μ kat/L	0.770	0.615	—	0.925	0.052
BS-200E ³	U/L	45.7	36.4	—	55.0	3.1	BS-480 ¹²	U/L	45.1	36.1	—	54.1	3.0
	μ kat/L	0.763	0.608	—	0.919	0.052		μ kat/L	0.753	0.603	—	0.903	0.050
BS-230 ⁴	U/L	46.7	37.4	—	56.0	3.1	BS-600 ¹³	U/L	45.0	36.0	—	54.0	3.0
	μ kat/L	0.780	0.625	—	0.935	0.052		μ kat/L	0.752	0.601	—	0.902	0.050
BS-240E ⁵	U/L	45.3	36.3	—	54.3	3.0	BS-600M ¹⁴	U/L	45.6	36.3	—	54.9	3.1
	μ kat/L	0.757	0.606	—	0.907	0.050		μ kat/L	0.762	0.606	—	0.917	0.052
BS-300 ⁶	U/L	44.5	35.5	—	53.5	3.0	BS-800 ¹⁵	U/L	45.0	36.0	—	54.0	3.0
	μ kat/L	0.743	0.593	—	0.893	0.050		μ kat/L	0.752	0.601	—	0.902	0.050
BS-330E ⁷	U/L	45.7	36.4	—	55.0	3.1	BS-2000 ¹⁶	U/L	45.2	36.2	—	54.2	3.0
	μ kat/L	0.763	0.608	—	0.919	0.052		μ kat/L	0.755	0.605	—	0.905	0.050
BS-360E ⁸	U/L	44.9	35.9	—	53.9	3.0	BS-2800M ¹⁷	U/L	45.2	36.2	—	54.2	3.0
	μ kat/L	0.750	0.600	—	0.900	0.050		μ kat/L	0.755	0.605	—	0.905	0.050
BS-380 ⁹	U/L	45.4	36.4	—	54.4	3.0							
	μ kat/L	0.758	0.608	—	0.908	0.050							

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

LIP	
English	Lipase
Русский	Липаза
Português	Lipase
Español	Lipasa
Italiano	Lipasi
Türkçe	Lipaz