

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray  
 Português: Planilha de valores do Sistema de Medição BS da Mindray  
 Español: Hoja de valores del sistema de medición Mindray BS  
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray  
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-300: BS-300, BS-320;
- 7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 9.BS-380: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

- 10.BS-400: BS-400, BS-420;
- 11.BS-430: BS-430, BS-450,BS-460;
- 12.BS-480: BS-480, BS-490;
- 13.BS-600: BS-600, BS-620;
- 14.BS-600M: BS-600M, BS-620M;
- 15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 17.BS-2800M:BS-2600M, BS-2800M;
- 18.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

LOT: 044824004

Expires: 2025-07-31

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

## LIP

English	Lipase
Русский	Липаза
Português	Lipase
Español	Lipasa
Italiano	Lipasi
Türkçe	Lipaz

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
LIP	U/L	BS-120 <sup>1</sup>	43.9	35.2 — 52.6	2.9	BS-400 <sup>10</sup>	45.2	36.2 — 54.2	3.0
		BS-200 <sup>2</sup>	43.8	35.1 — 52.5	2.9	BS-430 <sup>11</sup>	43.5	34.8 — 52.2	2.9
		BS-200E <sup>3</sup>	45.1	36.1 — 54.1	3.0	BS-480 <sup>12</sup>	44.0	35.3 — 52.7	2.9
		BS-240 <sup>4</sup>	44.9	35.9 — 53.9	3.0	BS-600 <sup>13</sup>	44.0	35.3 — 52.7	2.9
		BS-240E <sup>5</sup>	44.0	35.3 — 52.7	2.9	BS-600M <sup>14</sup>	44.2	35.2 — 53.2	3.0
		BS-300 <sup>6</sup>	42.9	34.2 — 51.6	2.9	BS-800 <sup>15</sup>	44.0	35.3 — 52.7	2.9
		BS-330E <sup>7</sup>	45.1	36.1 — 54.1	3.0	BS-2000 <sup>16</sup>	43.8	35.1 — 52.5	2.9
		BS-360E <sup>8</sup>	44.3	35.3 — 53.3	3.0	BS-2800M <sup>17</sup>	44.3	35.3 — 53.3	3.0
		BS-380 <sup>9</sup>	45.2	36.2 — 54.2	3.0				
				BS-120 <sup>1</sup>	0.733	0.588 — 0.878	0.048	BS-400 <sup>10</sup>	0.755
µkat/L		BS-200 <sup>2</sup>	0.731	0.586 — 0.877	0.048	BS-430 <sup>11</sup>	0.726	0.581 — 0.872	0.048
		BS-200E <sup>3</sup>	0.753	0.603 — 0.903	0.050	BS-480 <sup>12</sup>	0.735	0.590 — 0.880	0.048
		BS-240 <sup>4</sup>	0.750	0.600 — 0.900	0.050	BS-600 <sup>13</sup>	0.735	0.590 — 0.880	0.048
		BS-240E <sup>5</sup>	0.735	0.590 — 0.880	0.048	BS-600M <sup>14</sup>	0.738	0.588 — 0.888	0.050

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value $\pm$ 3SD)		1 SD	Model	Assay Value	Range(Assay Value $\pm$ 3SD)		1 SD		
		<b>BS-300</b> <sup>6</sup>	0.716	0.571	—	0.862	0.048	<b>BS-800</b> <sup>15</sup>	0.735	0.590	—	0.880	0.048
		<b>BS-330E</b> <sup>7</sup>	0.753	0.603	—	0.903	0.050	<b>BS-2000</b> <sup>16</sup>	0.731	0.586	—	0.877	0.048
		<b>BS-360E</b> <sup>8</sup>	0.740	0.590	—	0.890	0.050	<b>BS-2800M</b> <sup>17</sup>	0.740	0.590	—	0.890	0.050
		<b>BS-380</b> <sup>9</sup>	0.755	0.605	—	0.905	0.050						