

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-230**: BS-230, BS-240;

5. **BS-400**: BS-400, BS-420;

2. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-430**: BS-430, BS-450, BS-460;

3. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

7. **BS-480**: BS-480, BS-490;

BS-330E(V35.00) (Serial Number starts with "W8-" and

8. **BS-600**: BS-600, BS-620;

software

9. **BS-600M**: BS-600M, BS-620M;

version starts with "35.00"), **BS-350E(V35.00)**(Serial Number

10. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

starts with "W9-"and software version starts with "35.00");

BS-1800, BS-1800plus;

4. **BS-380**: BS-380, BS-390;

11. **BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

12. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 053524001

EXP: 2025-04-16

Abbreviated name	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range(Assay Value±3SD)		1 SD		
UIBC	BS-230 ¹	μmol/L	32.3	25.7	—	38.9	2.2	BS-480 ⁷	μmol/L	34.8	27.9	—	41.7	2.3
		μg/dL	181	144	—	217	12		μg/dL	195	156	—	233	13
	BS-240E ²	μmol/L	34.7	27.8	—	41.6	2.3	BS-600 ⁸	μmol/L	32.8	26.2	—	39.4	2.2
		μg/dL	194	155	—	233	13		μg/dL	183	146	—	220	12
	BS-360E ³	μmol/L	32.5	25.9	—	39.1	2.2	BS-600M ⁹	μmol/L	32.7	26.1	—	39.3	2.2
		μg/dL	182	145	—	219	12		μg/dL	183	146	—	220	12
	BS-380 ⁴	μmol/L	33.5	26.9	—	40.1	2.2	BS-800 ¹⁰	μmol/L	32.8	26.2	—	39.4	2.2
		μg/dL	187	150	—	224	12		μg/dL	183	146	—	220	12
	BS-400 ⁵	μmol/L	33.5	26.9	—	40.1	2.2	BS-2000 ¹¹	μmol/L	30.7	24.4	—	37.0	2.1
		μg/dL	187	150	—	224	12		μg/dL	172	136	—	207	12
	BS-430 ⁶	μmol/L	32.8	26.2	—	39.4	2.2	BS-2800M ¹²	μmol/L	31.1	24.8	—	37.4	2.1
		μg/dL	183	146	—	220	12		μg/dL	174	139	—	209	12

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон (результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa (Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore diconcentrazione ±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

UIBC

English	Unsaturated Iron Binding Capacity
Русский	ненасыщенная железосвязывающая способность
Português	Capacidade de ligação de ferro insaturado
Español	Capacidad de unión de hierro no saturado
Italiano	Capacità di legame del ferro insaturo
Türkçe	Unsature Demir Bağlama Kapasitesi