

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

6. **BS-400**: BS-400,BS-420;

2. **BS-230**: BS-230,BS-240;

7. **BS-430**: BS-430,BS-450,BS-460;

3. **BS-240E**: BS-240E, BS-240Pro;

8. **BS-480**: BS-480,BS-490;

4. **BS-360E**: BS-360E, BS-370E,BS-350S,BS-360S,

9. **BS-600**: BS-600, BS-620;

BS-330E (V35.00)Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5. **BS-380**: BS-380, BS-390;

12. **BS-2000**: BS-2000,BS-2200,BS-2000M,BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 054623007

EXP : 2024-12-18

Abbreviated name	RBP Control					1 SD
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		
RBP	BS-200E ¹	mg/dL	4.68	3.75	– 5.61	0.31
	BS-230 ²	mg/dL	4.68	3.75	– 5.61	0.31
	BS-240E ³	mg/dL	4.74	3.78	– 5.70	0.32
	BS-360E ⁴	mg/dL	4.74	3.78	– 5.70	0.32
	BS-380 ⁵	mg/dL	4.68	3.75	– 5.61	0.31
	BS-400 ⁶	mg/dL	4.68	3.75	– 5.61	0.31
	BS-430 ⁷	mg/dL	4.74	3.78	– 5.70	0.32
	BS-480 ⁸	mg/dL	4.74	3.78	– 5.70	0.32
	BS-600 ⁹	mg/dL	4.74	3.78	– 5.70	0.32
	BS-600M ¹⁰	mg/dL	4.74	3.78	– 5.70	0.32
	BS-800 ¹¹	mg/dL	4.74	3.78	– 5.70	0.32
	BS-2000 ¹²	mg/dL	4.85	3.89	– 5.81	0.32
	BS-2800M ¹³	mg/dL	4.74	3.78	– 5.70	0.32

English	Abbreviated name	RBP Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль РСБ	модель	Прибор	Результат анализа	Диапазон (результат анализа \pm 3CO)
Português	Nome abreviado	Controle RBP	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control RBP	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo RBP	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	RBP Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

RBP	
English	Retinol Binding Protein
Русский	Ретинол-связывающий белок
Português	Proteína transportadora do retinol
Español	Proteína fijadora de retinol
Italiano	Proteina legante il retinolo
Türkçe	Retinol Bağlayıcı Protein