

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS-240E, BS-240Pro;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

(Serial Number starts with "W9-" and software version starts with "35.00");

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

7. **BS-380**: BS-380, BS-390;

8. **BS-400**: BS-400, BS-420;

9. **BS-480**: BS-480, BS-490;

10. **BS-600**: BS-600, BS-620;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 057522001

: 2023-01-18

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	BS-120 ¹	U/L	1027	821	—	1233	69	BS-120 ¹	U/L	2225	1778	—	2672	149
		µkat/L	17.2	13.7	—	20.6	1.2		µkat/L	37.2	29.7	—	44.6	2.5
	BS-200 ²	U/L	1044	834	—	1254	70	BS-200 ²	U/L	2222	1775	—	2669	149
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	37.1	29.6	—	44.6	2.5
	BS-200E ³	U/L	867	693	—	1041	58	BS-200E ³	U/L	1899	1517	—	2281	127
		µkat/L	14.5	11.6	—	17.4	1.0		µkat/L	31.7	25.3	—	38.1	2.1
	BS-230 ⁴	U/L	975	779	—	1171	65	BS-230 ⁴	U/L	2115	1690	—	2540	142
		µkat/L	16.3	13.0	—	19.6	1.1		µkat/L	35.3	28.2	—	42.4	2.4
	BS-240E ⁵	U/L	971	776	—	1166	65	BS-240E ⁵	U/L	2117	1691	—	2543	142
		µkat/L	16.2	13.0	—	19.5	1.1		µkat/L	35.4	28.2	—	42.5	2.4
	BS-360E ⁶	U/L	944	754	—	1134	63	BS-360E ⁶	U/L	1928	1540	—	2316	129
		µkat/L	15.8	12.6	—	18.9	1.1		µkat/L	32.2	25.7	—	38.7	2.2
	BS-380 ⁷	U/L	963	769	—	1157	65	BS-380 ⁷	U/L	2079	1661	—	2497	139
		µkat/L	16.1	12.8	—	19.3	1.1		µkat/L	34.7	27.7	—	41.7	2.3
	BS-400 ⁸	U/L	941	752	—	1130	63	BS-400 ⁸	U/L	2059	1645	—	2473	138
		µkat/L	15.7	12.6	—	18.9	1.1		µkat/L	34.4	27.5	—	41.3	2.3
	BS-480 ⁹	U/L	987	789	—	1185	66	BS-480 ⁹	U/L	2095	1674	—	2516	140
		µkat/L	16.5	13.2	—	19.8	1.1		µkat/L	35.0	28.0	—	42.0	2.3

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	987	789	—	1185	66	BS-600 ¹⁰	U/L	2095	1674	—	2516	140
		μ kat/L	16.5	13.2	—	19.8	1.1		μ kat/L	35.0	28.0	—	42.0	2.3
	BS-800 ¹¹	U/L	987	789	—	1185	66	BS-800 ¹¹	U/L	2095	1674	—	2516	140
		μ kat/L	16.5	13.2	—	19.8	1.1		μ kat/L	35.0	28.0	—	42.0	2.3
	BS-2000 ¹²	U/L	1022	817	—	1227	68	BS-2000 ¹²	U/L	2142	1711	—	2573	144
		μ kat/L	17.1	13.6	—	20.5	1.1		μ kat/L	35.8	28.6	—	43.0	2.4
	BS-2800M ¹³	U/L	1028	821	—	1235	69	BS-2800M ¹³	U/L	2169	1733	—	2605	145
		μ kat/L	17.2	13.7	—	20.6	1.2		μ kat/L	36.2	28.9	—	43.5	2.4

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz