

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

(Serial Number starts with "W9-" and software version starts with "35.00");

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057522003

: 2023-09-29

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)			1 SD
G6PD	BS-120 ¹	U/L	1043	833	—	1253	70	BS-120 ¹	U/L	2263	1808	—	2718	152
		µkat/L	17.4	13.9	—	20.9	1.2		µkat/L	37.8	30.2	—	45.4	2.5
	BS-200 ²	U/L	1071	856	—	1286	72	BS-200 ²	U/L	2282	1823	—	2741	153
		µkat/L	17.9	14.3	—	21.5	1.2		µkat/L	38.1	30.4	—	45.8	2.6
	BS-200E ³	U/L	885	707	—	1063	59	BS-200E ³	U/L	1974	1577	—	2371	132
		µkat/L	14.8	11.8	—	17.8	1.0		µkat/L	33.0	26.3	—	39.6	2.2
	BS-230 ⁴	U/L	967	773	—	1161	65	BS-230 ⁴	U/L	2069	1653	—	2485	139
		µkat/L	16.1	12.9	—	19.4	1.1		µkat/L	34.6	27.6	—	41.5	2.3
	BS-240E ⁵	U/L	1032	825	—	1239	69	BS-240E ⁵	U/L	2175	1738	—	2612	146
		µkat/L	17.2	13.8	—	20.7	1.2		µkat/L	36.3	29.0	—	43.6	2.4
	BS-360E ⁶	U/L	937	749	—	1125	63	BS-360E ⁶	U/L	2013	1608	—	2418	135
		µkat/L	15.6	12.5	—	18.8	1.1		µkat/L	33.6	26.9	—	40.4	2.3
	BS-380 ⁷	U/L	977	781	—	1173	65	BS-380 ⁷	U/L	2111	1687	—	2535	141
		µkat/L	16.3	13.0	—	19.6	1.1		µkat/L	35.3	28.2	—	42.3	2.4
	BS-400 ⁸	U/L	934	746	—	1122	63	BS-400 ⁸	U/L	2056	1643	—	2469	138
		µkat/L	15.6	12.5	—	18.7	1.1		µkat/L	34.3	27.4	—	41.2	2.3
	BS-480 ⁹	U/L	1031	824	—	1238	69	BS-480 ⁹	U/L	2153	1720	—	2586	144
		µkat/L	17.2	13.8	—	20.7	1.2		µkat/L	36.0	28.7	—	43.2	2.4

Abbreviated name	G6PD Control (L)							G6PD Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
G6PD	BS-600 ¹⁰	U/L	1031	824	—	1238	69	BS-600 ¹⁰	U/L	2153	1720	—	2586	144
		μ kat/L	17.2	13.8	—	20.7	1.2		μ kat/L	36.0	28.7	—	43.2	2.4
	BS-600M ¹¹	U/L	/	/	—	/	/	BS-600M ¹¹	U/L	/	/	—	/	/
		μ kat/L	/	/	—	/	/		μ kat/L	/	/	—	/	/
	BS-800 ¹²	U/L	1031	824	—	1238	69	BS-800 ¹²	U/L	2153	1720	—	2586	144
		μ kat/L	17.2	13.8	—	20.7	1.2		μ kat/L	36.0	28.7	—	43.2	2.4
	BS-2000 ¹³	U/L	1060	847	—	1273	71	BS-2000 ¹³	U/L	2192	1751	—	2633	147
		μ kat/L	17.7	14.1	—	21.3	1.2		μ kat/L	36.6	29.2	—	44.0	2.5
	BS-2800M ¹⁴	U/L	1047	837	—	1257	70	BS-2800M ¹⁴	U/L	2186	1747	—	2625	146
		μ kat/L	17.5	14.0	—	21.0	1.2		μ kat/L	36.5	29.2	—	43.8	2.4

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz