

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

BS-1800, BS-1800plus;

(Serial Number starts with "W9-" and software version

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

starts with "35.00");

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057523001

: 2024-05-25

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	BS-120 ¹	U/L	1137	908	—	1366	76	BS-120 ¹	U/L	2419	1933	—	2905	162
		µkat/L	19.0	15.2	—	22.8	1.3		µkat/L	40.4	32.3	—	48.5	2.7
	BS-200 ²	U/L	1133	905	—	1361	76	BS-200 ²	U/L	2402	1919	—	2885	161
		µkat/L	18.9	15.1	—	22.7	1.3		µkat/L	40.1	32.0	—	48.2	2.7
	BS-200E ³	U/L	972	777	—	1167	65	BS-200E ³	U/L	2064	1649	—	2479	138
		µkat/L	16.2	13.0	—	19.5	1.1		µkat/L	34.5	27.5	—	41.4	2.3
	BS-230 ⁴	U/L	1010	807	—	1213	68	BS-230 ⁴	U/L	2120	1694	—	2546	142
		µkat/L	16.9	13.5	—	20.3	1.1		µkat/L	35.4	28.3	—	42.5	2.4
	BS-240E ⁵	U/L	1046	836	—	1256	70	BS-240E ⁵	U/L	2232	1783	—	2681	150
		µkat/L	17.5	14.0	—	21.0	1.2		µkat/L	37.3	29.8	—	44.8	2.5
	BS-360E ⁶	U/L	984	786	—	1182	66	BS-360E ⁶	U/L	2072	1656	—	2488	139
		µkat/L	16.4	13.1	—	19.7	1.1		µkat/L	34.6	27.7	—	41.5	2.3
	BS-380 ⁷	U/L	1002	801	—	1203	67	BS-380 ⁷	U/L	2153	1720	—	2586	144
		µkat/L	16.7	13.4	—	20.1	1.1		µkat/L	36.0	28.7	—	43.2	2.4
	BS-400 ⁸	U/L	1010	807	—	1213	68	BS-400 ⁸	U/L	2128	1700	—	2556	143
		µkat/L	16.9	13.5	—	20.3	1.1		µkat/L	35.5	28.4	—	42.7	2.4
BS-480 ⁹	U/L	1049	838	—	1260	70	BS-480 ⁹	U/L	2189	1749	—	2629	147	
	µkat/L	17.5	14.0	—	21.0	1.2		µkat/L	36.6	29.2	—	43.9	2.5	

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
G6PD	BS-600 ¹⁰	U/L	1049	838	—	1260	70	BS-600 ¹⁰	U/L	2189	1749	—	2629	147
		μ kat/L	17.5	14.0	—	21.0	1.2		μ kat/L	36.6	29.2	—	43.9	2.5
	BS-600M ¹¹	U/L	1036	828	—	1244	69	BS-600M ¹¹	U/L	2217	1771	—	2663	149
		μ kat/L	17.3	13.8	—	20.8	1.2		μ kat/L	37.0	29.6	—	44.5	2.5
	BS-800 ¹²	U/L	1049	838	—	1260	70	BS-800 ¹²	U/L	2189	1749	—	2629	147
		μ kat/L	17.5	14.0	—	21.0	1.2		μ kat/L	36.6	29.2	—	43.9	2.5
	BS-2000 ¹³	U/L	1089	870	—	1308	73	BS-2000 ¹³	U/L	2258	1804	—	2712	151
		μ kat/L	18.2	14.5	—	21.8	1.2		μ kat/L	37.7	30.1	—	45.3	2.5
	BS-2800M ¹⁴	U/L	1081	864	—	1298	72	BS-2800M ¹⁴	U/L	2241	1791	—	2691	150
		μ kat/L	18.1	14.4	—	21.7	1.2		μ kat/L	37.4	29.9	—	44.9	2.5

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz