

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-480**: BS-480, BS-490;

4. **BS-230**: BS-230, BS-240;

10. **BS-600**: BS-600, BS-620;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600M**: BS-600M, BS-620M;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

12. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)**

BS-1800, BS-1800plus;

(Serial Number starts with "W9-" and software version starts with "35.00");

13. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

14. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057523002

: 2024-08-22

Abbreviated name	G6PD Control (L)						G6PD Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
G6PD	BS-120 ¹	U/L	1077	861	—	1293	72	BS-120 ¹	U/L	2290	1830	—	2750	153
		µkat/L	18.0	14.4	—	21.6	1.2		µkat/L	38.2	30.6	—	45.9	2.6
	BS-200 ²	U/L	1107	884	—	1330	74	BS-200 ²	U/L	2346	1874	—	2818	157
		µkat/L	18.5	14.8	—	22.2	1.2		µkat/L	39.2	31.3	—	47.1	2.6
	BS-200E ³	U/L	956	764	—	1148	64	BS-200E ³	U/L	2014	1609	—	2419	135
		µkat/L	16.0	12.8	—	19.2	1.1		µkat/L	33.6	26.9	—	40.4	2.3
	BS-230 ⁴	U/L	1023	817	—	1229	69	BS-230 ⁴	U/L	2151	1719	—	2583	144
		µkat/L	17.1	13.6	—	20.5	1.2		µkat/L	35.9	28.7	—	43.1	2.4
	BS-240E ⁵	U/L	1068	853	—	1283	72	BS-240E ⁵	U/L	2238	1788	—	2688	150
		µkat/L	17.8	14.2	—	21.4	1.2		µkat/L	37.4	29.9	—	44.9	2.5
	BS-360E ⁶	U/L	1001	800	—	1202	67	BS-360E ⁶	U/L	2092	1672	—	2512	140
		µkat/L	16.7	13.4	—	20.1	1.1		µkat/L	34.9	27.9	—	42.0	2.3
	BS-380 ⁷	U/L	1025	819	—	1231	69	BS-380 ⁷	U/L	2157	1723	—	2591	145
		µkat/L	17.1	13.7	—	20.6	1.2		µkat/L	36.0	28.8	—	43.3	2.4
	BS-400 ⁸	U/L	1012	809	—	1215	68	BS-400 ⁸	U/L	2137	1707	—	2567	143
		µkat/L	16.9	13.5	—	20.3	1.1		µkat/L	35.7	28.5	—	42.9	2.4
	BS-480 ⁹	U/L	1055	843	—	1267	71	BS-480 ⁹	U/L	2174	1737	—	2611	146
		µkat/L	17.6	14.1	—	21.2	1.2		µkat/L	36.3	29.0	—	43.6	2.4

Abbreviated name	G6PD Control (L)							G6PD Control (H)						
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)			1 SD
G6PD	BS-600 ¹⁰	U/L	1055	843	—	1267	71	BS-600 ¹⁰	U/L	2174	1737	—	2611	146
		μ kat/L	17.6	14.1	—	21.2	1.2		μ kat/L	36.3	29.0	—	43.6	2.4
	BS-600M ¹¹	U/L	1066	852	—	1280	71	BS-600M ¹¹	U/L	2220	1774	—	2666	149
		μ kat/L	17.8	14.2	—	21.4	1.2		μ kat/L	37.1	29.6	—	44.5	2.5
	BS-800 ¹²	U/L	1055	843	—	1267	71	BS-800 ¹²	U/L	2174	1737	—	2611	146
		μ kat/L	17.6	14.1	—	21.2	1.2		μ kat/L	36.3	29.0	—	43.6	2.4
	BS-2000 ¹³	U/L	1094	874	—	1314	73	BS-2000 ¹³	U/L	2227	1779	—	2675	149
		μ kat/L	18.3	14.6	—	21.9	1.2		μ kat/L	37.2	29.7	—	44.7	2.5
	BS-2800M ¹⁴	U/L	1094	874	—	1314	73	BS-2800M ¹⁴	U/L	2249	1797	—	2701	151
		μ kat/L	18.3	14.6	—	21.9	1.2		μ kat/L	37.6	30.0	—	45.1	2.5

English	Abbreviated name	G6PD Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль G6PD	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle G6PD	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control G6PD	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo G6PD	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	G6PD Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

G6PD	
English	Glucose-6-Phosphate Dehydrogenase
Русский	глюкозо-6-фосфатдегидрогеназа
Português	Glicose-6-fosfato desidrogenase
Español	Glucosa-6-fosfato deshidrogenasa
Italiano	Glucosio-6-fosfato deidrogenasi
Türkçe	Glikoz 6 Fosfat Dehidrogenaz