

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов измерительной системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадает во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

7. **BS-380**: BS-380, BS-390;

2. **BS-200**: BS-200, BS-220;

8. **BS-400**: BS-400, BS-420;

3. **BS-200E**: BS-200E, BS-220E;

9. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-230**: BS-230, BS-240;

10. **BS-480**: BS-480, BS-490;

5. **BS-240E**: BS-240E, BS-240Pro;

11. **BS-600**: BS-600, BS-620;

6. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**

12. **BS-600M**: BS-600M, BS-620M;

(Serial Number starts with "W8-" and software version starts

13. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

with "35.00"), **BS-350E(V35.00)**(Serial Number starts with

BS-1800, BS-1800plus;

"W9-" and software version starts with "35.00");

14. **BS-2000**: BS-2000, BS-2200, BS-2000M,

BS-2200M;

15. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 057723009

EXP : 2024-11-13

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
ACE	BS-120 ¹	U/L	55.2	41.4	—	69.0	4.6	BS-120 ¹	U/L	158	119	—	197	13
		µkat/L	0.922	0.691	—	1.152	0.077		µkat/L	2.64	1.99	—	3.29	0.22
	BS-200 ²	U/L	52.6	39.4	—	65.8	4.4	BS-200 ²	U/L	158	119	—	197	13
		µkat/L	0.878	0.658	—	1.099	0.073		µkat/L	2.64	1.99	—	3.29	0.22
	BS-200E ³	U/L	50.1	37.5	—	62.7	4.2	BS-200E ³	U/L	145	109	—	181	12
		µkat/L	0.837	0.626	—	1.047	0.070		µkat/L	2.42	1.82	—	3.02	0.20
	BS-230 ⁴	U/L	55.2	41.4	—	69.0	4.6	BS-230 ⁴	U/L	158	119	—	197	13
		µkat/L	0.922	0.691	—	1.152	0.077		µkat/L	2.64	1.99	—	3.29	0.22
	BS-240E ⁵	U/L	44.5	33.4	—	55.6	3.7	BS-240E ⁵	U/L	134	101	—	167	11
		µkat/L	0.743	0.558	—	0.929	0.062		µkat/L	2.24	1.69	—	2.79	0.18
	BS-360E ⁶	U/L	44.5	33.4	—	55.6	3.7	BS-360E ⁶	U/L	134	101	—	167	11
		µkat/L	0.743	0.558	—	0.929	0.062		µkat/L	2.24	1.69	—	2.79	0.18
BS-380 ⁷	U/L	45.5	34.1	—	56.9	3.8	BS-380 ⁷	U/L	134	101	—	167	11	
	µkat/L	0.760	0.569	—	0.950	0.063		µkat/L	2.24	1.69	—	2.79	0.18	
BS-400 ⁸	U/L	46.3	34.9	—	57.7	3.8	BS-400 ⁸	U/L	133	100	—	166	11	
	µkat/L	0.773	0.583	—	0.964	0.063		µkat/L	2.22	1.67	—	2.77	0.18	

Abbreviated name	ACE Control (L)						ACE Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ACE	BS-430 ⁹	U/L	44.5	33.4	—	55.6	3.7	BS-430 ⁹	U/L	134	101	—	167	11
		μ kat/L	0.743	0.558	—	0.929	0.062		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-480 ¹⁰	U/L	44.5	33.4	—	55.6	3.7	BS-480 ¹⁰	U/L	134	101	—	167	11
		μ kat/L	0.743	0.558	—	0.929	0.062		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-600 ¹¹	U/L	44.5	33.4	—	55.6	3.7	BS-600 ¹¹	U/L	134	101	—	167	11
		μ kat/L	0.743	0.558	—	0.929	0.062		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-600M ¹²	U/L	46.7	35.0	—	58.4	3.9	BS-600M ¹²	U/L	141	105	—	177	12
		μ kat/L	0.780	0.585	—	0.975	0.065		μ kat/L	2.35	1.75	—	2.96	0.20
	BS-800 ¹³	U/L	44.5	33.4	—	55.6	3.7	BS-800 ¹³	U/L	134	101	—	167	11
		μ kat/L	0.743	0.558	—	0.929	0.062		μ kat/L	2.24	1.69	—	2.79	0.18
	BS-2000 ¹⁴	U/L	45.7	34.3	—	57.1	3.8	BS-2000 ¹⁴	U/L	136	103	—	169	11
		μ kat/L	0.763	0.573	—	0.954	0.063		μ kat/L	2.27	1.72	—	2.82	0.18
BS-2800M ¹⁵	U/L	48.2	36.2	—	60.2	4.0	BS-2800M ¹⁵	U/L	142	106	—	178	12	
	μ kat/L	0.805	0.605	—	1.005	0.067		μ kat/L	2.37	1.77	—	2.97	0.20	

English	Abbreviated name	ACE Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль ACE	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle ACE	Modelo	Unidade	Valores da análise	Faixa (Valores da análise \pm 3SD)
Español	nombre abreviado	Control ACE	modelo	Unidad	Valor de ensayo	Rango (Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo ACE	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	ACE Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ACE	
English	Angiotensin Converting Enzyme
Русский	ангиотензинпревращающего фермента
Português	Enzima de conversão da angiotensina
Español	Enzima convertidora de angiotensina
Italiano	Enzima convertitore dell'angiotensina
Türkçe	Anjiotensin Dönüştürücü Enzim