

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-480: BS-480, BS-490;

2.BS-240: BS-230, BS-240;

9.BS-600: BS-600, BS-620;

3.BS-240E: BS240E, BS240Pro;

10.BS-600M: BS-600M, BS-620M;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-"

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

and software version starts with "35.00");

13.BS-2800M:BS-2600M, BS-2800M;

5.BS-380: BS-380, BS-390;

14.For applicable analyzers of the analyte, please subject to the

6.BS-400: BS-400, BS-420;

parameter sheet and instrument.

7.BS-430: BS-430, BS-450,BS-460;

LOT : 058024002

🕒: 2025-06-11

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kisaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

FER			IgE		
English		Ferritin			Immunoglobulin E
Русский		ферритин			Иммуноглобулин E
Português		Ferritina			Imunoglobulina E
Español		Ferritina			Inmunoglobulina E
Italiano		Ferritina			Immunoglobulina E
Türkçe		Ferritin			İmmünglobulin E

MYO		
English		Myoglobin
Русский		Миоглобин
Português		Mioglobina
Español		Mioglobina
Italiano		Mioglobina
Türkçe		Myoglobin

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ng/mL		BS-200E ¹	124	106	—	142	6	BS-200E ¹	466	397	—	535	23
		BS-240 ²	124	106	—	142	6	BS-240 ²	475	403	—	547	24
		BS-240E ³	124	106	—	142	6	BS-240E ³	468	399	—	537	23
		BS-360E ⁴	124	106	—	142	6	BS-360E ⁴	468	399	—	537	23
		BS-380 ⁵	124	106	—	142	6	BS-380 ⁵	466	397	—	535	23
		BS-400 ⁶	124	106	—	142	6	BS-400 ⁶	466	397	—	535	23
		BS-430 ⁷	124	106	—	142	6	BS-430 ⁷	468	399	—	537	23
		BS-480 ⁸	124	106	—	142	6	BS-480 ⁸	468	399	—	537	23
		BS-600 ⁹	124	106	—	142	6	BS-600 ⁹	468	399	—	537	23
		BS-600M ¹⁰	123	105	—	141	6	BS-600M ¹⁰	457	388	—	526	23
		BS-800 ¹¹	124	106	—	142	6	BS-800 ¹¹	468	399	—	537	23

FER	pmol/L	BS-2000 ¹²	121	103	—	139	6	BS-2000 ¹²	451	382	—	520	23
		BS-2800M ¹³	123	105	—	141	6	BS-2800M ¹³	457	388	—	526	23
		BS-200E ¹	279	238	—	319	13	BS-200E ¹	1047	892	—	1202	52
		BS-240 ²	279	238	—	319	13	BS-240 ²	1067	906	—	1229	54
		BS-240E ³	279	238	—	319	13	BS-240E ³	1052	897	—	1207	52
		BS-360E ⁴	279	238	—	319	13	BS-360E ⁴	1052	897	—	1207	52
		BS-380 ⁵	279	238	—	319	13	BS-380 ⁵	1047	892	—	1202	52
		BS-400 ⁶	279	238	—	319	13	BS-400 ⁶	1047	892	—	1202	52
		BS-430 ⁷	279	238	—	319	13	BS-430 ⁷	1052	897	—	1207	52
		BS-480 ⁸	279	238	—	319	13	BS-480 ⁸	1052	897	—	1207	52
		BS-600 ⁹	279	238	—	319	13	BS-600 ⁹	1052	897	—	1207	52
		BS-600M ¹⁰	276	236	—	317	13	BS-600M ¹⁰	1027	872	—	1182	52
		BS-800 ¹¹	279	238	—	319	13	BS-800 ¹¹	1052	897	—	1207	52
		BS-2000 ¹²	272	231	—	312	13	BS-2000 ¹²	1013	858	—	1168	52
BS-2800M ¹³	276	236	—	317	13	BS-2800M ¹³	1027	872	—	1182	52		

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
IgE	IU/mL	BS-240 ²	79.2	55.5 — 102.9	7.9	BS-240 ²	529	370 — 688	53
		BS-240E ³	81.8	57.2 — 106.4	8.2	BS-240E ³	521	365 — 677	52
		BS-360E ⁴	82.6	57.7 — 107.5	8.3	BS-360E ⁴	507	354 — 660	51
		BS-380 ⁵	77.4	54.3 — 100.5	7.7	BS-380 ⁵	521	365 — 677	52
		BS-400 ⁶	80.3	56.3 — 104.3	8.0	BS-400 ⁶	527	368 — 686	53
		BS-480 ⁸	81.6	57.0 — 106.2	8.2	BS-480 ⁸	523	367 — 679	52
		BS-600 ⁹	81.6	57.0 — 106.2	8.2	BS-600 ⁹	523	367 — 679	52
		BS-600M ¹⁰	82.5	57.6 — 107.4	8.3	BS-600M ¹⁰	523	367 — 679	52
		BS-800 ¹¹	81.6	57.0 — 106.2	8.2	BS-800 ¹¹	523	367 — 679	52
		BS-2000 ¹²	78.5	54.8 — 102.2	7.9	BS-2000 ¹²	517	361 — 673	52
		BS-2800M ¹³	82.5	57.6 — 107.4	8.3	BS-2800M ¹³	523	367 — 679	52

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
MYO	ng/mL	BS-240 ²	105	90 — 120	5	BS-240 ²	403	343 — 463	20
		BS-240E ³	106	91 — 121	5	BS-240E ³	409	349 — 469	20
		BS-360E ⁴	106	91 — 121	5	BS-360E ⁴	409	349 — 469	20
		BS-380 ⁵	106	91 — 121	5	BS-380 ⁵	407	347 — 467	20
		BS-400 ⁶	106	91 — 121	5	BS-400 ⁶	407	347 — 467	20
		BS-430 ⁷	106	91 — 121	5	BS-430 ⁷	399	339 — 459	20
		BS-480 ⁸	106	91 — 121	5	BS-480 ⁸	409	349 — 469	20
		BS-600 ⁹	106	91 — 121	5	BS-600 ⁹	409	349 — 469	20
		BS-600M ¹⁰	106	91 — 121	5	BS-600M ¹⁰	407	347 — 467	20
		BS-800 ¹¹	106	91 — 121	5	BS-800 ¹¹	409	349 — 469	20
		BS-2000 ¹²	108	93 — 123	5	BS-2000 ¹²	420	357 — 483	21
		BS-2800M ¹³	106	91 — 121	5	BS-2800M ¹³	407	347 — 467	20

mindray

Multimmun Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 058024002

🕒 2025-06-11

