

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

5. **BS-380**: BS-380, BS-390;

2. **BS-230**: BS-230, BS-240;

6. **BS-400**: BS-400, BS-420;

3. **BS-240E**: BS240E, BS240Pro;

7. **BS-430**: BS-430, BS-450, BS-460;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

8. **BS-480**: BS-480, BS-490;

BS-330E(V35.00)(Serial Number starts with "W8"-

9. **BS-600**: BS-600, BS-620;

and software version starts with "35.00"),

10. **BS-600M**: BS-600M, BS-620M;

BS-350E(V35.00) (Serial Number starts with "W9"-

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

and software version starts with "35.00");

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 059022007

EXP: 2024-02-29

Abbreviated name	Cystatin C Control (L)					Cystatin C Control (H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
CysC II	BS-200E ¹	mg/L	0.683	0.581 — 0.785	0.034	BS-200E ¹	mg/L	1.15	0.98 — 1.32	0.06
	BS-230 ²	mg/L	0.677	0.575 — 0.779	0.034	BS-230 ²	mg/L	1.15	0.98 — 1.32	0.06
	BS-240E ³	mg/L	0.670	0.570 — 0.771	0.034	BS-240E ³	mg/L	1.14	0.97 — 1.31	0.06
	BS-360E ⁴	mg/L	0.684	0.581 — 0.787	0.034	BS-360E ⁴	mg/L	1.15	0.98 — 1.32	0.06
	BS-380 ⁵	mg/L	0.683	0.581 — 0.785	0.034	BS-380 ⁵	mg/L	1.15	0.98 — 1.32	0.06
	BS-400 ⁶	mg/L	0.676	0.575 — 0.777	0.034	BS-400 ⁶	mg/L	1.14	0.97 — 1.31	0.06
	BS-430 ⁷	mg/L	0.684	0.581 — 0.787	0.034	BS-430 ⁷	mg/L	1.15	0.98 — 1.32	0.06
	BS-480 ⁸	mg/L	0.675	0.574 — 0.776	0.034	BS-480 ⁸	mg/L	1.13	0.96 — 1.30	0.06
	BS-600 ⁹	mg/L	0.684	0.581 — 0.787	0.034	BS-600 ⁹	mg/L	1.15	0.98 — 1.32	0.06
	BS-600M ¹⁰	mg/L	0.701	0.596 — 0.806	0.035	BS-600M ¹⁰	mg/L	1.16	0.99 — 1.33	0.06
	BS-800 ¹¹	mg/L	0.684	0.581 — 0.787	0.034	BS-800 ¹¹	mg/L	1.15	0.98 — 1.32	0.06
	BS-2000 ¹²	mg/L	0.690	0.587 — 0.794	0.035	BS-2000 ¹²	mg/L	1.17	0.99 — 1.35	0.06
	BS-2800M ¹³	mg/L	0.699	0.594 — 0.804	0.035	BS-2800M ¹³	mg/L	1.15	0.98 — 1.32	0.06

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль цистатина С	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

CysC	
English	Cystatin C
Русский	цистатина С
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C