

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

5. **BS-380**: BS-380, BS-390;

2. **BS-230**: BS-230, BS-240;

6. **BS-400**: BS-400, BS-420;

3. **BS-240E**: BS240E, BS240Pro;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

9. **BS-600**: BS-600, BS-620;

BS-330E(V35.00)(Serial Number starts with "W8-"

10. **BS-600M**: BS-600M, BS-620M;

and software version starts with "35.00"),

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

BS-350E(V35.00) (Serial Number starts with "W9-"

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

and software version starts with "35.00");

13. **BS-2800M**: BS-2600M, BS-2800M.

LOT : 059023005

EXP : 2024-09-29

Abbreviated name	Cystatin C Control (L)					Cystatin C Control (H)				
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)	1 SD
CysC II	BS-200E ¹	mg/L	0.742	0.631 — 0.853	0.037	BS-200E ¹	mg/L	1.29	1.10 — 1.48	0.06
	BS-230 ²	mg/L	0.746	0.634 — 0.858	0.037	BS-230 ²	mg/L	1.29	1.10 — 1.48	0.06
	BS-240E ³	mg/L	0.750	0.638 — 0.863	0.038	BS-240E ³	mg/L	1.29	1.10 — 1.48	0.06
	BS-360E ⁴	mg/L	0.754	0.641 — 0.867	0.038	BS-360E ⁴	mg/L	1.31	1.11 — 1.51	0.07
	BS-380 ⁵	mg/L	0.742	0.631 — 0.853	0.037	BS-380 ⁵	mg/L	1.29	1.10 — 1.48	0.06
	BS-400 ⁶	mg/L	0.747	0.635 — 0.859	0.037	BS-400 ⁶	mg/L	1.28	1.09 — 1.47	0.06
	BS-430 ⁷	mg/L	0.754	0.641 — 0.867	0.038	BS-430 ⁷	mg/L	1.31	1.11 — 1.51	0.07
	BS-480 ⁸	mg/L	0.745	0.633 — 0.857	0.037	BS-480 ⁸	mg/L	1.29	1.10 — 1.48	0.06
	BS-600 ⁹	mg/L	0.754	0.641 — 0.867	0.038	BS-600 ⁹	mg/L	1.31	1.11 — 1.51	0.07
	BS-600M ¹⁰	mg/L	0.745	0.633 — 0.857	0.037	BS-600M ¹⁰	mg/L	1.30	1.11 — 1.50	0.07
	BS-800 ¹¹	mg/L	0.754	0.641 — 0.867	0.038	BS-800 ¹¹	mg/L	1.31	1.11 — 1.51	0.07
	BS-2000 ¹²	mg/L	0.771	0.655 — 0.887	0.039	BS-2000 ¹²	mg/L	1.31	1.11 — 1.51	0.07
	BS-2800M ¹³	mg/L	0.749	0.637 — 0.861	0.037	BS-2800M ¹³	mg/L	1.29	1.10 — 1.48	0.06

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

CysC	
English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C