

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray

Português : Planilha de valores do Sistema de Medição BS da Mindray

Español : Hoja de valores del sistema de medición Mindray BS

Italiano : Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

1.**BS-200E**: BS-200E, BS-220E;

2.**BS-230**: BS-230, BS-240;

3.**BS-240E**: BS240E, BS240Pro;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-"

and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with"W9-"

and software version starts with "35.00");

5.**BS-380**: BS-380, BS-390;

6.**BS-400**: BS-400, BS-420;

7.**BS-430**: BS-430, BS-450, BS-460;

8.**BS-480**: BS-480, BS-490;

9.**BS-600**: BS-600, BS-620;

10.**BS-600M**: BS-600M, BS-620M;

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

LOT : 059023009

EXP : 2025-03-12

Abbreviated name	Cystatin C Control (L)						Cystatin C Control (H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
CysC II	BS-200E ¹	mg/L	0.721	0.613	—	0.829	0.036	BS-200E ¹	mg/L	1.23	1.05	—	1.41	0.06
	BS-230 ²	mg/L	0.691	0.586	—	0.796	0.035	BS-230 ²	mg/L	1.24	1.06	—	1.42	0.06
	BS-240E ³	mg/L	0.761	0.647	—	0.875	0.038	BS-240E ³	mg/L	1.28	1.10	—	1.46	0.06
	BS-360E ⁴	mg/L	0.726	0.618	—	0.834	0.036	BS-360E ⁴	mg/L	1.25	1.07	—	1.43	0.06
	BS-380 ⁵	mg/L	0.721	0.613	—	0.829	0.036	BS-380 ⁵	mg/L	1.23	1.05	—	1.41	0.06
	BS-400 ⁶	mg/L	0.712	0.604	—	0.820	0.036	BS-400 ⁶	mg/L	1.24	1.06	—	1.42	0.06
	BS-430 ⁷	mg/L	0.726	0.618	—	0.834	0.036	BS-430 ⁷	mg/L	1.25	1.07	—	1.43	0.06
	BS-480 ⁸	mg/L	0.724	0.616	—	0.832	0.036	BS-480 ⁸	mg/L	1.24	1.06	—	1.42	0.06
	BS-600 ⁹	mg/L	0.726	0.618	—	0.834	0.036	BS-600 ⁹	mg/L	1.25	1.07	—	1.43	0.06
	BS-600M ¹⁰	mg/L	0.755	0.641	—	0.869	0.038	BS-600M ¹⁰	mg/L	1.24	1.06	—	1.42	0.06
	BS-800 ¹¹	mg/L	0.726	0.618	—	0.834	0.036	BS-800 ¹¹	mg/L	1.25	1.07	—	1.43	0.06
	BS-2000 ¹²	mg/L	0.730	0.619	—	0.841	0.037	BS-2000 ¹²	mg/L	1.26	1.08	—	1.44	0.06
	BS-2800M ¹³	mg/L	0.750	0.636	—	0.864	0.038	BS-2800M ¹³	mg/L	1.24	1.06	—	1.42	0.06

English	Abbreviated name	Cystatin C Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль цистатина C	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controlo Cistatina C	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control cistatina C	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo Cistatina C	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	Cystatin C Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

CysC	
English	Cystatin C
Русский	цистатина C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C