

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-480: BS-480, BS-490;

2.BS-240: BS-230, BS-240;

9.BS-600: BS-600, BS-620;

3.BS-240E: BS240E, BS240Pro;

10.BS-600M: BS-600M, BS-620M;

4.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

5.BS-380: BS-380, BS-390;

13.BS-2800M:BS-2600M, BS-2800M;

6.BS-400: BS-400, BS-420;

14.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

7.BS-430: BS-430, BS-450,BS-460;

LOT : 059024003

: 2025-07-05

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

CysC

English	Cystatin C
Русский	цистатин C
Português	Cistatina C
Español	cistatina C
Italiano	cistatina C
Türkçe	Cystatin C

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
CysC II	mg/L	BS-200E ¹	0.727	0.619	—	0.835	0.036	BS-200E ¹	1.21	1.03	—	1.39	0.06
		BS-240 ²	0.728	0.620	—	0.836	0.036	BS-240 ²	1.21	1.03	—	1.39	0.06
		BS-240E ³	0.778	0.661	—	0.895	0.039	BS-240E ³	1.26	1.08	—	1.44	0.06
		BS-360E ⁴	0.742	0.631	—	0.853	0.037	BS-360E ⁴	1.22	1.04	—	1.40	0.06
		BS-380 ⁵	0.727	0.619	—	0.835	0.036	BS-380 ⁵	1.21	1.03	—	1.39	0.06
		BS-400 ⁶	0.732	0.621	—	0.843	0.037	BS-400 ⁶	1.22	1.04	—	1.40	0.06
		BS-430 ⁷	0.742	0.631	—	0.853	0.037	BS-430 ⁷	1.22	1.04	—	1.40	0.06
		BS-480 ⁸	0.731	0.620	—	0.842	0.037	BS-480 ⁸	1.21	1.03	—	1.39	0.06
		BS-600 ⁹	0.742	0.631	—	0.853	0.037	BS-600 ⁹	1.22	1.04	—	1.40	0.06
		BS-600M ¹⁰	0.777	0.660	—	0.894	0.039	BS-600M ¹⁰	1.23	1.05	—	1.41	0.06
		BS-800 ¹¹	0.742	0.631	—	0.853	0.037	BS-800 ¹¹	1.22	1.04	—	1.40	0.06
		BS-2000 ¹²	0.753	0.639	—	0.867	0.038	BS-2000 ¹²	1.24	1.06	—	1.42	0.06
		BS-2800M ¹³	0.756	0.642	—	0.870	0.038	BS-2800M ¹³	1.21	1.03	—	1.39	0.06

mindray
Cystatin C Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 059024003

EXP 2025-07-05

