

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
 Português : Planilha de valores do Sistema de Medição BS da Mindray
 Español : Hoja de valores del sistema de medición Mindray BS
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский : Данные совпадают во всех группах.

Português : A dados de cada grupo é a mesma.

Español : la datos de cada grupo es la misma.

Italiano : la dati di ogni gruppo è la stessa.

Türkçe : her grubun veri aynıdır.

1.**BS-200E**: BS-200E, BS-220E;

6.**BS-400**: BS-400, BS-420;

2.**BS-230**: BS-230, BS-240;

7.**BS-430**: BS-430, BS-450, BS-460;

3.**BS-240E**: BS-240E, BS-240Pro;

8.**BS-480**: BS-480, BS-490;

4.**BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

9.**BS-600**: BS-600, BS-620;

BS-330E(V35.00)(Serial Number starts with "W8-"
and software version starts with "35.00");

10.**BS-600M**: BS-600M, BS-620M;

BS-350E(V35.00) (Serial Number starts with "W9-"
and software version starts with "35.00");

11.**BS-800**: BS-800, BS-820, BS-800M, BS-820M,
BS-1800, BS-1800plus;

5.**BS-380**: BS-380, BS-390;

12.**BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13.**BS-2800M**:BS-2600M, BS-2800M.

14.**Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.**

LOT : 059222006

EXP : 2024-04-10

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	109	71	—	147	13	BS-480 ⁸	109	71	—	147	13
		BS-230 ²	109	71	—	147	13	BS-600 ⁹	109	71	—	147	13
		BS-240E ³	109	71	—	147	13	BS-600M ¹⁰	103	67	—	139	12
		BS-360E ⁴	109	71	—	147	13	BS-800 ¹¹	109	71	—	147	13
		BS-380 ⁵	109	71	—	147	13	BS-2000 ¹²	105	68	—	142	12
		BS-400 ⁶	109	71	—	147	13	BS-2800M ¹³	103	67	—	139	12
		BS-430 ⁷	109	71	—	147	13						
RF II	IU/mL	BS-200E ¹	29.8	19.3	—	40.3	3.5	BS-480 ⁸	29.6	19.2	—	40.0	3.5
		BS-230 ²	30.5	19.8	—	41.2	3.6	BS-600 ⁹	30.1	19.5	—	40.7	3.5
		BS-240E ³	29.1	18.9	—	39.3	3.4	BS-600M ¹⁰	29.4	19.1	—	39.7	3.4
		BS-360E ⁴	29.6	19.2	—	40.0	3.5	BS-800 ¹¹	29.6	19.2	—	40.0	3.5
		BS-380 ⁵	29.9	19.4	—	40.4	3.5	BS-2000 ¹²	29.5	19.1	—	39.9	3.5
		BS-400 ⁶	29.9	19.4	—	40.4	3.5	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	29.6	19.2	—	40.0	3.5						
RF III	IU/mL	BS-800 ¹¹	40.5	26.3	—	54.7	4.7	BS-2800M ¹³	40.0	26.0	—	54.0	4.7
		BS-2000 ¹²	41.2	26.7	—	55.7	4.8						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray
 Português: Planilha de valores do Sistema de Medição BS da Mindray
 Español: Hoja de valores del sistema de medición Mindray BS
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

1. **BS-200E**: BS-200E, BS-220E;

2. **BS-230**: BS-230, BS-240;

3. **BS-240E**: BS-240E, BS-240Pro;

4. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S,

BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"),

BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

5. **BS-380**: BS-380, BS-390;

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

Español: la datos de cada grupo es la misma.

6. **BS-400**: BS-400, BS-420;

7. **BS-430**: BS-430, BS-450, BS-460;

8. **BS-480**: BS-480, BS-490;

9. **BS-600**: BS-600, BS-620;

10. **BS-600M**: BS-600M, BS-620M;

11. **BS-800**: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

12. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

13. **BS-2800M**: BS-2600M, BS-2800M.

14. Control target value of RF reagent kit corresponds to RF III, and control target value of RF II reagent kit corresponds to RF II.

LOT: 059222006

: 2024-04-10

Abbreviated name	Unit	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ASO II	IU/mL	BS-200E ¹	291	189	—	393	34	BS-480 ⁸	286	186	—	386	33
		BS-230 ²	283	184	—	382	33	BS-600 ⁹	286	186	—	386	33
		BS-240E ³	286	186	—	386	33	BS-600M ¹⁰	280	182	—	378	33
		BS-360E ⁴	286	186	—	386	33	BS-800 ¹¹	286	186	—	386	33
		BS-380 ⁵	291	189	—	393	34	BS-2000 ¹²	281	182	—	380	33
		BS-400 ⁶	291	189	—	393	34	BS-2800M ¹³	280	182	—	378	33
		BS-430 ⁷	286	186	—	386	33						
RF II	IU/mL	BS-200E ¹	78.6	51.0	—	106.2	9.2	BS-480 ⁸	72.5	47.1	—	97.9	8.5
		BS-230 ²	74.5	48.4	—	100.6	8.7	BS-600 ⁹	74.3	48.2	—	100.4	8.7
		BS-240E ³	76.6	49.7	—	103.5	9.0	BS-600M ¹⁰	71.7	46.5	—	96.9	8.4
		BS-360E ⁴	72.5	47.1	—	97.9	8.5	BS-800 ¹¹	72.5	47.1	—	97.9	8.5
		BS-380 ⁵	73.3	47.6	—	99.0	8.6	BS-2000 ¹²	74.0	48.0	—	100.0	8.7
		BS-400 ⁶	73.3	47.6	—	99.0	8.6	BS-2800M ¹³	/	/	—	/	/
		BS-430 ⁷	72.5	47.1	—	97.9	8.5						
RF III	IU/mL	BS-800 ¹¹	81.7	53.0	—	110.4	9.6	BS-2800M ¹³	81.7	53.0	—	110.4	9.6
		BS-2000 ¹²	82.3	53.4	—	111.2	9.6						

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ASO			RF		
English	Antistreptolysin "O"		Rheumatoid Factor		
Русский	антистрептолизина O		ревматоидные факторы		
Português	Antiestreptolisina "O"		Fator Reumatoide		
Español	antiestreptolisina "O"		factor reumatoide		
Italiano	Anti-Streptolisina "O"		fattore reumatoide		
Türkçe	Antistreptolisin "O"		Romatoid Faktör		