

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1. **BS-120**: BS-120, BS-130;

2. **BS-200**: BS-200, BS-220;

3. **BS-200E**: BS-200E, BS-220E;

4. **BS-230**: BS-230, BS-240;

5. **BS-240E**: BS240E, BS240Pro;

6. **BS-330**: BS-330, BS-350;

7. **BS-330E**: BS-330E (Serial Number starts with "XQ-"),

BS-350E (Serial Number starts with "XS-");

8. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)** (Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

9. **BS-380**: BS-380, BS-390;

10. **BS-400**: BS-400, BS-420;

11. **BS-430**: BS-430, BS-450, BS-460;

12. **BS-480**: BS-480, BS-490;

13. **BS-600**: BS-600, BS-620;

14. **BS-600M**: BS-600M, BS-620M;

15. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

16. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M;

17. **BS-2800M**: BS-2600M, BS-2800M.

LOT: 060522005

EXP: 2023-11-27

Abbreviated name	TPUC Control (L)					TPUC Control (H)								
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
TPUC	BS-120 ¹	mg/L	127	108	—	146	6	BS-120 ¹	mg/L	481	409	—	553	24
	BS-200 ²	mg/L	127	108	—	146	6	BS-200 ²	mg/L	481	409	—	553	24
	BS-200E ³	mg/L	126	107	—	145	6	BS-200E ³	mg/L	494	420	—	568	25
	BS-230 ⁴	mg/L	127	108	—	146	6	BS-230 ⁴	mg/L	475	404	—	546	24
	BS-240E ⁵	mg/L	129	110	—	148	6	BS-240E ⁵	mg/L	498	423	—	573	25
	BS-330 ⁶	mg/L	127	108	—	146	6	BS-330 ⁶	mg/L	481	409	—	553	24
	BS-330E ⁷	mg/L	126	107	—	145	6	BS-330E ⁷	mg/L	494	420	—	568	25
	BS-360E ⁸	mg/L	129	110	—	148	6	BS-360E ⁸	mg/L	498	423	—	573	25
	BS-380 ⁹	mg/L	121	103	—	139	6	BS-380 ⁹	mg/L	485	412	—	558	24
	BS-400 ¹⁰	mg/L	123	105	—	141	6	BS-400 ¹⁰	mg/L	486	413	—	559	24
	BS-430 ¹¹	mg/L	129	110	—	148	6	BS-430 ¹¹	mg/L	498	423	—	573	25
	BS-480 ¹²	mg/L	129	110	—	148	6	BS-480 ¹²	mg/L	498	423	—	573	25
	BS-600 ¹³	mg/L	125	106	—	144	6	BS-600 ¹³	mg/L	492	418	—	566	25
	BS-600M ¹⁴	mg/L	129	110	—	148	6	BS-600M ¹⁴	mg/L	494	420	—	568	25
	BS-800 ¹⁵	mg/L	129	110	—	148	6	BS-800 ¹⁵	mg/L	498	423	—	573	25
	BS-2000 ¹⁶	mg/L	132	112	—	152	7	BS-2000 ¹⁶	mg/L	511	434	—	588	26
	BS-2800M ¹⁷	mg/L	129	110	—	148	6	BS-2800M ¹⁷	mg/L	494	420	—	568	25

English	Abbreviated name	TPUC Control	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	Контроль общего белка в моче/СМЖ	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Controle PTUL	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	Control TPUC	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	Controllo TPUC	modelli	Unità	Valori di dosaggio	Intervallo(valore di concentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	TPUC Kontrolü	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

TPUC	
English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)