

**Value sheet of Mindray BS Measurement System**

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-430: BS-430, BS-450,BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

12.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

13.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

15.BS-1000M: BS-1000M, BS-1100M;

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

17.BS-2800M:BS-2600M, BS-2800M;

8.BS-380: BS-380, BS-390;

18.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

9.BS-400: BS-400, BS-420;

**LOT** : 060524003

**EXP**: 2025-07-11

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**TPUC**

English	Total Protein in Urine/CSF (TPUC)
Русский	Общий белок в моче/СМЖ (TPUC)
Português	Proteína total em urina e LCR (PTUL)
Español	proteína total en orina/LCR (TPUC)
Italiano	Proteine Totali Urinarie/CSF (TPUC)
Türkçe	BOS/İdrar Protein (TPUC)

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
TPUC	mg/L	BS-120 <sup>1</sup>	120	102	—	138	6	BS-120 <sup>1</sup>	461	392	—	530	23
		BS-200 <sup>2</sup>	120	102	—	138	6	BS-200 <sup>2</sup>	461	392	—	530	23
		BS-200E <sup>3</sup>	125	107	—	143	6	BS-200E <sup>3</sup>	497	422	—	572	25
		BS-240 <sup>4</sup>	120	102	—	138	6	BS-240 <sup>4</sup>	461	392	—	530	23
		BS-240E <sup>5</sup>	125	107	—	143	6	BS-240E <sup>5</sup>	497	422	—	572	25
		BS-330E <sup>6</sup>	125	107	—	143	6	BS-330E <sup>6</sup>	497	422	—	572	25
		BS-360E <sup>7</sup>	125	107	—	143	6	BS-360E <sup>7</sup>	497	422	—	572	25
		BS-380 <sup>8</sup>	120	102	—	138	6	BS-380 <sup>8</sup>	495	420	—	570	25
		BS-400 <sup>9</sup>	125	107	—	143	6	BS-400 <sup>9</sup>	495	420	—	570	25
		BS-430 <sup>10</sup>	125	107	—	143	6	BS-430 <sup>10</sup>	497	422	—	572	25
		BS-480 <sup>11</sup>	125	107	—	143	6	BS-480 <sup>11</sup>	495	420	—	570	25
		BS-600 <sup>12</sup>	125	107	—	143	6	BS-600 <sup>12</sup>	495	420	—	570	25
		BS-600M <sup>13</sup>	125	107	—	143	6	BS-600M <sup>13</sup>	495	420	—	570	25
		BS-800 <sup>14</sup>	125	107	—	143	6	BS-800 <sup>14</sup>	497	422	—	572	25
		BS-1000M <sup>15</sup>	125	107	—	143	6	BS-1000M <sup>15</sup>	495	420	—	570	25

---

<b>BS-2000<sup>16</sup></b>	128	110	—	146	6	<b>BS-2000<sup>16</sup></b>	502	427	—	577	25
<b>BS-2800M<sup>17</sup></b>	125	107	—	143	6	<b>BS-2800M<sup>17</sup></b>	495	420	—	570	25

---

# mindray

## TPUC Control

For use on: BS-1000M, BS-2000,  
BS-2800M, BS-600M,  
BS-800

**LOT** 060524003

**2025-07-11**

