

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

- | | |
|-------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. BS-120 : BS-120, BS-130; | 10. BS-360E : BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00) (Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00"); |
| 2. BS-180 : BS-180, BS-190; | 11. BS-380 : BS-380, BS-390; |
| 3. BS-200 : BS-200, BS-220; | 12. BS-400 : BS-400, BS-420; |
| 4. BS-200E : BS-200E, BS-220E; | 13. BS-430 : BS-410, BS-430, BS-450, BS-460, BS-470; |
| 5. BS-230 : BS-230, BS-240, BS-280; | 14. BS-480 : BS-480, BS-490; |
| 6. BS-240E : BS-240E, BS-240Pro; | 15. BS-600 : BS-600, BS-620; |
| 7. BS-300 : BS-300, BS-320; | 16. BS-600M : BS-600M, BS-620M; |
| 8. BS-330 : BS-330, BS-350; | 17. BS-800 : BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus; |
| 9. BS-330E : BS-330E (Serial Number starts with "XQ-"), BS-350E (Serial Number starts with "XS-"); | 18. BS-830 : BS-830, BS-830S, BS-840, BS-850, BS-860, BS-870; |
| | 19. BS-2000 : BS-2000, BS-2200, BS-2000M, BS-2200M. |
| | 20. BS-2800M : BS-2600M, BS-2800M |

LOT: 061522009

EXP: 2024-01-31

Abbreviated name	ADA Control(L)						ADA Control(H)							
	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)		1 SD		
ADA	BS-120	¹ U/L	27.4	22.5	—	32.3	1.6	BS-120	¹ U/L	123	101	—	145	7
	BS-180	² U/L	27.4	22.5	—	32.3	1.6	BS-180	² U/L	123	101	—	145	7
	BS-200	³ U/L	27.4	22.5	—	32.3	1.6	BS-200	³ U/L	127	104	—	150	8
	BS-200E	⁴ U/L	27.3	22.4	—	32.2	1.6	BS-200E	⁴ U/L	128	105	—	151	8
	BS-230	⁵ U/L	28.1	23.0	—	33.2	1.7	BS-230	⁵ U/L	127	104	—	150	8
	BS-240E	⁶ U/L	27.0	22.1	—	31.9	1.6	BS-240E	⁶ U/L	126	103	—	149	8
	BS-300	⁷ U/L	28.3	23.2	—	33.4	1.7	BS-300	⁷ U/L	129	106	—	152	8
	BS-330	⁸ U/L	27.4	22.5	—	32.3	1.6	BS-330	⁸ U/L	127	104	—	150	8
	BS-330E	⁹ U/L	27.3	22.4	—	32.2	1.6	BS-330E	⁹ U/L	128	105	—	151	8
	BS-360E	¹⁰ U/L	27.0	22.1	—	31.9	1.6	BS-360E	¹⁰ U/L	126	103	—	149	8
	BS-380	¹¹ U/L	27.2	22.3	—	32.1	1.6	BS-380	¹¹ U/L	127	104	—	150	8
	BS-400	¹² U/L	27.3	22.4	—	32.2	1.6	BS-400	¹² U/L	127	104	—	150	8
	BS-430	¹³ U/L	26.9	22.1	—	31.7	1.6	BS-430	¹³ U/L	126	103	—	149	8
	BS-480	¹⁴ U/L	27.0	22.1	—	31.9	1.6	BS-480	¹⁴ U/L	127	104	—	150	8
	BS-600	¹⁵ U/L	27.2	22.3	—	32.1	1.6	BS-600	¹⁵ U/L	127	104	—	150	8
	BS-600M	¹⁶ U/L	27.4	22.5	—	32.3	1.6	BS-600M	¹⁶ U/L	129	106	—	152	8
	BS-800	¹⁷ U/L	27.2	22.3	—	32.1	1.6	BS-800	¹⁷ U/L	127	104	—	150	8
	BS-830	¹⁸ U/L	27.7	22.7	—	32.7	1.7	BS-830	¹⁸ U/L	130	107	—	153	8
	BS-2000	¹⁹ U/L	27.5	22.6	—	32.5	1.7	BS-2000	¹⁹ U/L	129	106	—	152	8
	BS-2800M	²⁰ U/L	27.4	22.5	—	32.3	1.6	BS-2800M	²⁰ U/L	129	106	—	152	8

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value \pm 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа \pm 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise \pm 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo \pm 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo(valore diconcentrazione \pm 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri \pm 3SD)

ADA	
English	Adenosine Deaminase
Русский	является аденозиндезаминаза
Português	Adenosina deaminase
Español	adenosina desaminasa
Italiano	adenosina deamminasi
Türkçe	Adenozin Deaminazdır