

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1. **BS-120**: BS-120, BS-130;

10. **BS-360E**: BS-360E, BS-370E, BS-350S, BS-360S, **BS-330E(V35.00)**(Serial Number starts with "W8-" and software version starts with "35.00"), **BS-350E(V35.00)** (Serial Number starts with "W9-" and software version starts with "35.00");

2. **BS-180**: BS-180, BS-190;

11. **BS-380**: BS-380, BS-390;

3. **BS-200**: BS-200, BS-220;

12. **BS-400**: BS-400, BS-420;

4. **BS-200E**: BS-200E, BS-220E;

13. **BS-430**: BS-410, BS-430, BS-450, BS-460, BS-470;

5. **BS-230**: BS-230, BS-240, BS-280;

14. **BS-480**: BS-480, BS-490;

6. **BS-240E**: BS-240E, BS-240Pro;

15. **BS-600**: BS-600, BS-620;

7. **BS-300**: BS-300, BS-320;

16. **BS-600M**: BS-600M, BS-620M;

8. **BS-330**: BS-330, BS-350;

17. **BS-800**: BS-800, BS-820, BS-800M, BS-820M,

BS-1800, BS-1800plus;

18. **BS-830**: BS-830, BS-830S, BS-840, BS-850,

BS-860, BS-870;

9. **BS-330E**: **BS-330E**(Serial Number starts with "XQ-"),

19. **BS-2000**: BS-2000, BS-2200, BS-2000M, BS-2200M.

**BS-350E**(Serial Number starts with "XS-");

20. **BS-2800M**: BS-2600M, BS-2800M

LOT : 061523001

⌚ : 2024-06-26

Abbreviated name	ADA Control(L)						ADA Control(H)							
	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Unit	Assay Value	Range (Assay Value±3SD)		1 SD		
ADA	BS-120	<sup>1</sup> U/L	28.0	23.0	—	33.0	1.7	BS-120	<sup>1</sup> U/L	124	102	—	146	7
	BS-180	<sup>2</sup> U/L	28.0	23.0	—	33.0	1.7	BS-180	<sup>2</sup> U/L	124	102	—	146	7
	BS-200	<sup>3</sup> U/L	27.9	22.9	—	32.9	1.7	BS-200	<sup>3</sup> U/L	124	102	—	146	7
	BS-200E	<sup>4</sup> U/L	27.5	22.6	—	32.5	1.7	BS-200E	<sup>4</sup> U/L	126	103	—	149	8
	BS-230	<sup>5</sup> U/L	28.8	23.6	—	34.0	1.7	BS-230	<sup>5</sup> U/L	126	103	—	149	8
	BS-240E	<sup>6</sup> U/L	27.9	22.9	—	32.9	1.7	BS-240E	<sup>6</sup> U/L	125	103	—	148	8
	BS-300	<sup>7</sup> U/L	28.7	23.5	—	33.9	1.7	BS-300	<sup>7</sup> U/L	126	103	—	149	8
	BS-330	<sup>8</sup> U/L	27.9	22.9	—	32.9	1.7	BS-330	<sup>8</sup> U/L	124	102	—	146	7
	BS-330E	<sup>9</sup> U/L	27.5	22.6	—	32.5	1.7	BS-330E	<sup>9</sup> U/L	126	103	—	149	8
	BS-360E	<sup>10</sup> U/L	27.9	22.9	—	32.9	1.7	BS-360E	<sup>10</sup> U/L	125	103	—	148	8
	BS-380	<sup>11</sup> U/L	28.0	23.0	—	33.0	1.7	BS-380	<sup>11</sup> U/L	125	103	—	148	8
	BS-400	<sup>12</sup> U/L	27.8	22.8	—	32.8	1.7	BS-400	<sup>12</sup> U/L	125	103	—	148	8
	BS-430	<sup>13</sup> U/L	27.4	22.5	—	32.3	1.6	BS-430	<sup>13</sup> U/L	124	102	—	146	7
	BS-480	<sup>14</sup> U/L	27.8	22.8	—	32.8	1.7	BS-480	<sup>14</sup> U/L	126	103	—	149	8
	BS-600	<sup>15</sup> U/L	27.7	22.7	—	32.7	1.7	BS-600	<sup>15</sup> U/L	125	103	—	148	8
	BS-600M	<sup>16</sup> U/L	27.9	22.9	—	32.9	1.7	BS-600M	<sup>16</sup> U/L	127	104	—	150	8
	BS-800	<sup>17</sup> U/L	27.7	22.7	—	32.7	1.7	BS-800	<sup>17</sup> U/L	126	103	—	149	8
	BS-830	<sup>18</sup> U/L	28.1	23.0	—	33.2	1.7	BS-830	<sup>18</sup> U/L	127	104	—	150	8
	BS-2000	<sup>19</sup> U/L	28.2	23.1	—	33.3	1.7	BS-2000	<sup>19</sup> U/L	127	104	—	150	8
	BS-2800M	<sup>20</sup> U/L	27.9	22.9	—	32.9	1.7	BS-2800M	<sup>20</sup> U/L	127	104	—	150	8

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo(valore diconcentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

ADA	
English	Adenosine Deaminase
Русский	является аденозиндеаминаза
Português	Adenosina deaminase
Español	adenosina desaminasa
Italiano	adenosina deamminasi
Türkçe	Adenozin Deaminazdır