

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200: BS-200, BS-220, BS-330, BS-350;

8.BS-480: BS-480, BS-490;

2.BS-200E: BS-200E, BS-220E,;

9.BS-600: BS-600, BS-620;

3.BS-300: BS-300, BS-320;

10.BS-600M: BS-600M, BS-620M;

4.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

11.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

5.BS-380: BS-380, BS-390;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-400: BS-400, BS-420;

13.BS-2800M:BS-2600M, BS-2800M;

7.BS-430: BS-430, BS-450,BS-460;

14.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 061623008

🕒: 2025-02-01

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

AFU

English	α-L-Fucosidase
Русский	α-L-фукозидаза
Português	α-L-Fucosidase
Español	α-L-fucosidasa
Italiano	α-L-Fucosidasi
Türkçe	α-L-Fukosidazdır

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
AFU	U/L	BS-200 ¹	27.4	22.0	—	32.8	1.8	BS-200 ¹	72.0	57.6	—	86.4	4.8
		BS-200E ²	24.6	19.8	—	29.4	1.6	BS-200E ²	63.9	51.0	—	76.8	4.3
		BS-300 ³	26.0	20.9	—	31.1	1.7	BS-300 ³	67.9	54.4	—	81.4	4.5
		BS-330E ⁴	24.6	19.8	—	29.4	1.6	BS-330E ⁴	63.9	51.0	—	76.8	4.3
		BS-380 ⁵	25.1	20.0	—	30.2	1.7	BS-380 ⁵	66.4	53.2	—	79.6	4.4
		BS-400 ⁶	24.3	19.5	—	29.1	1.6	BS-400 ⁶	64.0	51.1	—	76.9	4.3
		BS-430 ⁷	26.8	21.4	—	32.2	1.8	BS-430 ⁷	69.0	55.2	—	82.8	4.6
		BS-480 ⁸	24.6	19.8	—	29.4	1.6	BS-480 ⁸	64.4	51.5	—	77.3	4.3
		BS-600 ⁹	24.6	19.8	—	29.4	1.6	BS-600 ⁹	63.7	50.8	—	76.6	4.3
		BS-600M ¹⁰	30.5	24.5	—	36.5	2.0	BS-600M ¹⁰	78.8	62.9	—	94.7	5.3
		BS-800 ¹¹	24.1	19.3	—	28.9	1.6	BS-800 ¹¹	62.6	50.0	—	75.2	4.2
		BS-2000 ¹²	27.6	22.2	—	33.0	1.8	BS-2000 ¹²	71.6	57.2	—	86.0	4.8
		BS-2800M ¹³	28.5	22.8	—	34.2	1.9	BS-2800M ¹³	74.1	59.1	—	89.1	5.0

mindray
AFU Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 061623008

2025-02-01

