

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

10.BS-400: BS-400, BS-420;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

11.BS-430: BS-430, BS-450,BS-460;

3.BS-200E: BS-200E, BS-220E,;

12.BS-480: BS-480, BS-490;

4.BS-240: BS-230, BS-240;

13.BS-600: BS-600, BS-620;

5.BS-240E: BS240E, BS240Pro;

14.BS-600M: BS-600M, BS-620M;

6.BS-300: BS-300, BS-320;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

17.BS-2800M:BS-2600M, BS-2800M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

18.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

9.BS-380: BS-380, BS-390;

LOT : 061824001

: 2025-07-25

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

CO2			TBA		
English	Carbon Dioxide	Total Bile Acids			
Русский	Двуокись углерода	желчные кислоты			
Português	Dióxido de Carbono	Ácido Biliar Total			
Español	Dióxido de carbono	ácidos biliares totales			
Italiano	Anidride carbonica	acidi biliari totali			
Türkçe	Karbon Dioksit	Total Safra Asitleridir			

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
CO2	mmol/L	BS-200 ²	12.8	10.1	—	15.5	0.9	BS-200 ²	19.3	15.4	—	23.2	1.3
		BS-200E ³	13.6	10.9	—	16.3	0.9	BS-200E ³	19.9	16.0	—	23.8	1.3
		BS-240 ⁴	13.0	10.3	—	15.7	0.9	BS-240 ⁴	18.9	15.0	—	22.8	1.3
		BS-240E ⁵	13.5	10.8	—	16.2	0.9	BS-240E ⁵	20.0	16.1	—	23.9	1.3
		BS-300 ⁶	13.1	10.4	—	15.8	0.9	BS-300 ⁶	19.6	15.7	—	23.5	1.3
		BS-330E ⁷	13.6	10.9	—	16.3	0.9	BS-330E ⁷	19.9	16.0	—	23.8	1.3
		BS-360E ⁸	13.5	10.8	—	16.2	0.9	BS-360E ⁸	20.0	16.1	—	23.9	1.3
		BS-380 ⁹	13.6	10.9	—	16.3	0.9	BS-380 ⁹	19.9	16.0	—	23.8	1.3
		BS-400 ¹⁰	13.6	10.9	—	16.3	0.9	BS-400 ¹⁰	19.8	15.9	—	23.7	1.3
		BS-430 ¹¹	13.5	10.8	—	16.2	0.9	BS-430 ¹¹	20.0	16.1	—	23.9	1.3
		BS-480 ¹²	13.5	10.8	—	16.2	0.9	BS-480 ¹²	20.0	16.1	—	23.9	1.3
		BS-600 ¹³	13.5	10.8	—	16.2	0.9	BS-600 ¹³	20.0	16.1	—	23.9	1.3
		BS-600M ¹⁴	14.0	11.3	—	16.7	0.9	BS-600M ¹⁴	20.0	16.1	—	23.9	1.3
		BS-800 ¹⁵	13.5	10.8	—	16.2	0.9	BS-800 ¹⁵	20.0	16.1	—	23.9	1.3
		BS-2000 ¹⁶	13.3	10.6	—	16.0	0.9	BS-2000 ¹⁶	20.0	16.1	—	23.9	1.3

BS-2800M¹⁷ 13.8 11.1 — 16.5 0.9 BS-2800M¹⁷ 20.3 16.1 — 24.5 1.4

Abbreviated name	Unit	Level: L/N				Level: H/P			
		Model	Assay Value	Range (Assay Value±3SD)	1 SD	Model	Assay Value	Range (Assay Value±3SD)	1 SD
TBA	μmol/L	BS-120 ¹	20.8	16.6 — 25.0	1.4	BS-120 ¹	34.4	27.5 — 41.3	2.3
		BS-200 ²	20.7	16.5 — 24.9	1.4	BS-200 ²	33.2	26.6 — 39.8	2.2
		BS-200E ³	21.3	17.1 — 25.5	1.4	BS-200E ³	33.7	26.8 — 40.6	2.3
		BS-240 ⁴	21.1	16.9 — 25.3	1.4	BS-240 ⁴	34.5	27.6 — 41.4	2.3
		BS-240E ⁵	21.2	17.0 — 25.4	1.4	BS-240E ⁵	34.4	27.5 — 41.3	2.3
		BS-300 ⁶	21.8	17.3 — 26.3	1.5	BS-300 ⁶	35.6	28.4 — 42.8	2.4
		BS-330E ⁷	21.3	17.1 — 25.5	1.4	BS-330E ⁷	33.7	26.8 — 40.6	2.3
		BS-360E ⁸	21.0	16.8 — 25.2	1.4	BS-360E ⁸	34.8	27.9 — 41.7	2.3
		BS-380 ⁹	21.7	17.2 — 26.2	1.5	BS-380 ⁹	34.8	27.9 — 41.7	2.3
		BS-400 ¹⁰	21.4	17.2 — 25.6	1.4	BS-400 ¹⁰	34.6	27.7 — 41.5	2.3
		BS-430 ¹¹	21.4	17.2 — 25.6	1.4	BS-430 ¹¹	34.4	27.5 — 41.3	2.3
		BS-480 ¹²	21.5	17.3 — 25.7	1.4	BS-480 ¹²	33.7	26.8 — 40.6	2.3
		BS-600 ¹³	21.2	17.0 — 25.4	1.4	BS-600 ¹³	34.5	27.6 — 41.4	2.3
		BS-600M ¹⁴	21.5	17.3 — 25.7	1.4	BS-600M ¹⁴	34.3	27.4 — 41.2	2.3
		BS-800 ¹⁵	21.8	17.3 — 26.3	1.5	BS-800 ¹⁵	35.5	28.3 — 42.7	2.4
		BS-2000 ¹⁶	21.9	17.4 — 26.4	1.5	BS-2000 ¹⁶	35.1	27.9 — 42.3	2.4
		BS-2800M ¹⁷	21.5	17.3 — 25.7	1.4	BS-2800M ¹⁷	34.3	27.4 — 41.2	2.3

mindray

CO2 and TBA Multi Control

For use on: BS-2000, BS-2800M,
BS-600M, BS-800

LOT 061824001

🕒 2025-07-25

