

## Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray  
 Português: Planilha de valores do Sistema de Medição BS da Mindray  
 Español: Hoja de valores del sistema de medición Mindray BS  
 Italiano: Scheda dei valori del sistema di misurazione BS di Mindray  
 Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português: A dados de cada grupo é a mesma.

Italiano: la dati di ogni gruppo è la stessa.

1.BS-200: BS-200, BS-220, BS-330, BS-350;

2.BS-200E: BS-200E, BS-220E,;

3.BS-240: BS-230, BS-240;

4.BS-240E: BS240E, BS240Pro;

5.BS-300: BS-300, BS-320;

6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

8.BS-380: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español: la datos de cada grupo es la misma.

Türkçe: her grubun veri aynıdır.

9.BS-400: BS-400, BS-420;

10.BS-430: BS-430, BS-450,BS-460;

11.BS-480: BS-480, BS-490;

12.BS-600: BS-600, BS-620;

13.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

14.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

15.BS-2800M:BS-2600M, BS-2800M;

16.For applicable analyzers of the analyte, please subject to the parameter sheet and instrument.

LOT: 145124006

2025-02-21

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

## CO2

English	Carbon Dioxide
Русский	Двуокись углерода
Português	Dióxido de Carbono
Español	Dióxido de carbono
Italiano	Anidride carbonica
Türkçe	Karbon Dioksit

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
mmol/L		BS-200 <sup>1</sup>	12.8	10.1 — 15.5	0.9	BS-400 <sup>9</sup>	13.6	10.9 — 16.3	0.9
		BS-200E <sup>2</sup>	13.6	10.9 — 16.3	0.9	BS-430 <sup>10</sup>	13.5	10.8 — 16.2	0.9
		BS-240 <sup>3</sup>	13.0	10.3 — 15.7	0.9	BS-480 <sup>11</sup>	13.5	10.8 — 16.2	0.9
		BS-240E <sup>4</sup>	13.5	10.8 — 16.2	0.9	BS-600 <sup>12</sup>	13.5	10.8 — 16.2	0.9
		BS-300 <sup>5</sup>	13.1	10.4 — 15.8	0.9	BS-800 <sup>13</sup>	13.5	10.8 — 16.2	0.9
		BS-330E <sup>6</sup>	13.6	10.9 — 16.3	0.9	BS-2000 <sup>14</sup>	13.3	10.6 — 16.0	0.9
		BS-360E <sup>7</sup>	13.5	10.8 — 16.2	0.9	BS-2800M <sup>15</sup>	13.8	11.1 — 16.5	0.9
CO2		BS-380 <sup>8</sup>	13.6	10.9 — 16.3	0.9				
		BS-200 <sup>1</sup>	0.0128	0.0101 — 0.0155	0.0009	BS-400 <sup>9</sup>	0.0136	0.0109 — 0.0163	0.0009
		BS-200E <sup>2</sup>	0.0136	0.0109 — 0.0163	0.0009	BS-430 <sup>10</sup>	0.0135	0.0108 — 0.0162	0.0009
		BS-240 <sup>3</sup>	0.0130	0.0103 — 0.0157	0.0009	BS-480 <sup>11</sup>	0.0135	0.0108 — 0.0162	0.0009
		BS-240E <sup>4</sup>	0.0135	0.0108 — 0.0162	0.0009	BS-600 <sup>12</sup>	0.0135	0.0108 — 0.0162	0.0009
		BS-300 <sup>5</sup>	0.0131	0.0104 — 0.0158	0.0009	BS-800 <sup>13</sup>	0.0135	0.0108 — 0.0162	0.0009
		BS-330E <sup>6</sup>	0.0136	0.0109 — 0.0163	0.0009	BS-2000 <sup>14</sup>	0.0133	0.0106 — 0.0160	0.0009
BS-360E <sup>7</sup>	0.0135	0.0108 — 0.0162	0.0009	BS-2800M <sup>15</sup>	0.0138	0.0111 — 0.0165	0.0009		

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
		<b>BS-380<sup>B</sup></b>	0.0136	0.0109 — 0.0163	0.0009				

# mindray

## CO2 Control N

For use on: BS-2000, BS-2800M,  
BS-800

LOT 145124006

2025-02-21

