

## Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray  
 Português : Planilha de valores do Sistema de Medição BS da Mindray  
 Español : Hoja de valores del sistema de medición Mindray BS  
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray  
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



**The data of each group is same.**

**Русский: Данные совпадают во всех группах.**

**Português : A dados de cada grupo é a mesma.**

**Español : la datos de cada grupo es la misma.**

**Italiano : la dati di ogni gruppo è la stessa.**

**Türkçe : her grubun veri aynıdır.**

- 1.BS-200: BS-200, BS-220, BS-330, BS-350;
- 2.BS-200E: BS-200E, BS-220E,;
- 3.BS-240: BS-230, BS-240;
- 4.BS-240E: BS240E, BS240Pro;
- 5.BS-300: BS-300, BS-320;
- 6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 8.BS-380: BS-380, BS-390;
- 9.BS-400: BS-400, BS-420;

- 10.BS-430: BS-430, BS-450,BS-460;
- 11.BS-480: BS-480, BS-490;
- 12.BS-600: BS-600, BS-620;
- 13.BS-600M: BS-600M, BS-620M;
- 14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 15.BS-1000M: BS-1000M, BS-1100M;
- 16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 17.BS-2800M:BS-2600M, BS-2800M;
- 18.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

**LOT : 145224004**

**🕒 : 2026-03-20**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
<b>Русский</b>	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
<b>Português</b>	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
<b>Español</b>	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
<b>Italiano</b>	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
<b>Türkçe</b>	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

**CHE**

<b>English</b>	Cholinesterase
<b>Русский</b>	Холинэстераза
<b>Português</b>	Colinesterase
<b>Español</b>	Colinesterasa
<b>Italiano</b>	Colinesterasi
<b>Türkçe</b>	Kolinesteraz

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
CHE	U/L	<b>BS-200</b> <sup>1</sup>	5884	4702 — 7066	394	<b>BS-430</b> <sup>10</sup>	5929	4738 — 7120	397
		<b>BS-200E</b> <sup>2</sup>	5887	4705 — 7069	394	<b>BS-480</b> <sup>11</sup>	5992	4789 — 7195	401
		<b>BS-240</b> <sup>3</sup>	5884	4702 — 7066	394	<b>BS-600</b> <sup>12</sup>	5992	4789 — 7195	401
		<b>BS-240E</b> <sup>4</sup>	5793	4629 — 6957	388	<b>BS-600M</b> <sup>13</sup>	5992	4789 — 7195	401
		<b>BS-300</b> <sup>5</sup>	5992	4789 — 7195	401	<b>BS-800</b> <sup>14</sup>	5929	4738 — 7120	397
	<b>BS-330E</b> <sup>6</sup>	5887	4705 — 7069	394	<b>BS-1000M</b> <sup>15</sup>	5992	4789 — 7195	401	
	<b>BS-360E</b> <sup>7</sup>	5887	4705 — 7069	394	<b>BS-2000</b> <sup>16</sup>	5960	4763 — 7157	399	
	<b>BS-380</b> <sup>8</sup>	5992	4789 — 7195	401	<b>BS-2800M</b> <sup>17</sup>	5992	4789 — 7195	401	
	<b>BS-400</b> <sup>9</sup>	5992	4789 — 7195	401					
	µkat/L	<b>BS-200</b> <sup>1</sup>	98.3	78.5 — 118.0	6.6	<b>BS-430</b> <sup>10</sup>	99.0	79.1 — 118.9	6.6
<b>BS-200E</b> <sup>2</sup>		98.3	78.6 — 118.1	6.6	<b>BS-480</b> <sup>11</sup>	100	80 — 120	7	
<b>BS-240</b> <sup>3</sup>		98.3	78.5 — 118.0	6.6	<b>BS-600</b> <sup>12</sup>	100	80 — 120	7	
<b>BS-240E</b> <sup>4</sup>		96.7	77.3 — 116.2	6.5	<b>BS-600M</b> <sup>13</sup>	100	80 — 120	7	
<b>BS-300</b> <sup>5</sup>		100	80 — 120	7	<b>BS-800</b> <sup>14</sup>	99.0	79.1 — 118.9	6.6	

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value $\pm$ 3SD)			1 SD	Model	Assay Value	Range(Assay Value $\pm$ 3SD)			1 SD
		<b>BS-330E<sup>6</sup></b>	98.3	78.6	—	118.1	6.6	<b>BS-1000M<sup>15</sup></b>	100	80	—	120	7
		<b>BS-360E<sup>7</sup></b>	98.3	78.6	—	118.1	6.6	<b>BS-2000<sup>16</sup></b>	99.5	79.5	—	119.5	6.7
		<b>BS-380<sup>8</sup></b>	100	80	—	120	7	<b>BS-2800M<sup>17</sup></b>	100	80	—	120	7
		<b>BS-400<sup>9</sup></b>	100	80	—	120	7						