

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
 Português : Planilha de valores do Sistema de Medição BS da Mindray
 Español : Hoja de valores del sistema de medición Mindray BS
 Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
 Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Português : A dados de cada grupo é a mesma.

Italiano : la dati di ogni gruppo è la stessa.

- 1.BS-120: BS-120, BS-130, BS-180, BS-190;
- 2.BS-200: BS-200, BS-220, BS-330, BS-350;
- 3.BS-200E: BS-200E, BS-220E,;
- 4.BS-240: BS-230, BS-240;
- 5.BS-240E: BS240E, BS240Pro;
- 6.BS-300: BS-300, BS-320;
- 7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")
- 8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");
- 9.BS-380: BS-380, BS-390;

Русский: Данные совпадают во всех группах.

Español : la datos de cada grupo es la misma.

Türkçe : her grubun veri aynıdır.

- 10.BS-400: BS-400, BS-420;
- 11.BS-430: BS-430, BS-450,BS-460;
- 12.BS-480: BS-480, BS-490;
- 13.BS-600: BS-600, BS-620;
- 14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;
- 15.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;
- 16.BS-2800M:BS-2600M, BS-2800M;
- 17.Applicable models of the chemistry shall be subject to the parameter list and instructions.

LOT : 145424001

🕒 : 2025-03-02

| English | Abbreviated name | Model | Unit | Assay Value | Range (Assay Value±3SD) |
|------------------|--------------------------|---------|---------|--------------------|--|
| Русский | сокращенное наименование | модель | Прибор | Результат анализа | Диапазон(результат анализа ± 3CO) |
| Português | Nome abreviado | Modelo | Unidade | Valores da análise | Faixa(Valores da análise ±3SD) |
| Español | nombre abreviado | modelo | Unidad | Valor de ensayo | Rango(Valor de ensayo ±3SD) |
| Italiano | abbreviazione | modelli | Unità | Valori di dosaggio | Intervallo (valore di concentrazione±3 DS) |
| Türkçe | kısaltılmış ad | model | Ünite | Tayin Değeri | Aralık (Tayin Değeri±3SD) |

FUN

| | |
|------------------|--------------|
| English | Fructosamine |
| Русский | Фруктозамин |
| Português | Frutosamina |
| Español | Fructosamina |
| Italiano | Fruttosamina |
| Türkçe | Fruktozamin |

| Abbreviated name | Unit | Model | Assay Value | Range(Assay Value±3SD) | 1 SD | Model | Assay Value | Range(Assay Value±3SD) | 1 SD |
|------------------|--------|----------------------------|-------------|------------------------|------|------------------------------|-------------|------------------------|------|
| FUN | µmol/L | BS-120¹ | 459 | 390 — 528 | 23 | BS-380⁹ | 444 | 378 — 510 | 22 |
| | | BS-200² | 455 | 386 — 524 | 23 | BS-400¹⁰ | 457 | 388 — 526 | 23 |
| | | BS-200E³ | 478 | 406 — 550 | 24 | BS-430¹¹ | 471 | 399 — 543 | 24 |
| | | BS-240⁴ | 455 | 386 — 524 | 23 | BS-480¹² | 467 | 398 — 536 | 23 |
| | | BS-240E⁵ | 464 | 395 — 533 | 23 | BS-600¹³ | 464 | 395 — 533 | 23 |
| | | BS-300⁶ | 467 | 398 — 536 | 23 | BS-800¹⁴ | 467 | 398 — 536 | 23 |
| | | BS-330E⁷ | 478 | 406 — 550 | 24 | BS-2000¹⁵ | 472 | 400 — 544 | 24 |
| | | BS-360E⁸ | 471 | 399 — 543 | 24 | BS-2800M¹⁶ | 456 | 387 — 525 | 23 |

mindray
FUN Control P

For use on: BS-2000, BS-2800M,
BS-800

LOT 145424001

2025-03-02

