

Coagulation Control

Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)		CX-9000, CX-9010, CX-9200 and CX-9210 Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	10.74	8.59 ~ 12.89
	8623050118		APTT	s	27.80	22.24 ~ 33.36
			TT	s	15.42	12.34 ~ 18.50
	2024-05-23		FIB	g/L	2.82	2.26 ~ 3.38
		mg/dL		282	226 ~ 338	

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation/Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)		C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	12.25	9.80 ~ 14.70
	8623050118		APTT	s	29.41	23.53 ~ 35.29
			TT	s	13.90	11.12 ~ 16.68
	2024-05-23		FIB	g/L	2.73	2.18 ~ 3.28
		mg/dL		273	218 ~ 328	