



BeneHeart D3 Defibrillator Accessories

| Power Cable | |
|--------------------------------|---|
| DA8K-10-14453 | UK Power Cord |
| ECG - Adult/Paediatric | |
| 0010-30-42719 | ECG Trunk Cable, 3/5- lead, Adu/Paed, 12 Pin, Defib-Proof, IEC (order lead wires separately) |
| 0010-30-42732 | ECG Lead Wires, 3- lead, Pinch, Adu/Paed, IEC, Long 1m |
| 0010-30-42733 | ECG Lead Wires, 3-lead, Snap, Adu/Paed, TPU, IEC, Long 1m |
| 0010-30-42730 | ECG Lead Wires, 5-lead, Pinch, Adu, TPU, IEC, Long 1m |
| 0010-30-42736 | ECG Lead Wires, 5-lead, Snap, Adu, TPU, IEC, Long 1m |
| ECG - Paediatric | |
| 0010-30-42720 | ECG Trunk Cable, 3-lead, Paed, 12 Pin, Defib-Proof, IEC, 3m (order lead wires separately) |
| 0010-30-42901 | ECG Lead Wires, 3-lead, Snap, Paed, TPU, IEC, 0.6m |
| 0010-30-42899 | ECG Lead Wires, 3-lead, Pinch, Paed, TPU, IEC, 0.6m |
| NIBP | |
| 6200-30-09688 | NIBP Hose, Adu/Paed/Inf, 3 m |
| 040-000978-00 | NIBP Cuff, Reusable, CM1302, Child, 18-26 cm, with Connector |
| 040-000983-00 | NIBP Cuff, Reusable, CM1303, Adult, 25-35 cm, with Connector |
| 040-000988-00 | NIBP Cuff, Reusable, CM1304, Large Adult, 33-47 cm, with Connector |
| 040-000968-00 | NIBP Cuff, Reusable, CM1300, Small Infant, 7-13 cm, with Connector |
| 040-000973-00 | NIBP Cuff, Reusable, CM1301, Infant, 10-19 cm, with Connector |
| 040-000993-00 | NIBP Cuff, Reusable, CM1305, Thigh 46-66cm, with Connector |
| Mindray SpO² | |
| 0010-20-42710 | Mindray SpO ² Interface Cable, 2.5m, 7 pin |
| 512F-30-28263 | Mindray 512F Reusable SpO ² Finger Sensor, Adult (all Mindray SpO ²) |
| 512H-30-79061 | Mindray 512H Reusable SpO ² Finger Sensor, Paediatric 10-30 kg (all Mindray SpO ²) |
| 009-005088-00 | 520P Single Patient Use Sensor, Paediatric 10kg to 50kg (20) |
| 512E-30-90390 | Mindray 512E Adult Reusable SpO ² Soft Finger Sensor |
| 512G-30-90607 | Mindray 512G Paediatric Reusable SpO ² Soft Finger Sensor |
| 115-034527-00 | SPO ² Sensor (Reusable, Adu/Paed, Ear-Clip) |

| | |
|---|---|
| Masimo SpO² | |
| 040-000332-00 | Masimo SpO ² Interface Cable, 8-pin connector |
| 115-004961-00 | Masimo SpO ² Accessory Kit (Adult), 8-pin connector |
| 040-005973-00 | Masimo RDSET Interface Cable - 8-pin connector |
| CO² | |
| 045-003134-00 | CO ² Microstream Adapter |
| 6800-30-50619 | Microstream CO ² kit, Adu |
| 040-000027-00 | Microstream CO ² kit, Paed |
| 040-000026-00 | Microstream CO ² kit, Neo |
| Defibrillator Electrodes (Pads) | |
| 115-001288-00 | Mindray multifunction electrode (Adult) |
| 115-035426-00 | Mindray multifunction electrode (Adult) Pack of 5 |
| 115-040517-00 | Mindray multifunction electrode Radio translucent (Adult) |
| 0651-30-77007 | Mindray multifunction electrode Radio translucent (Adult) Pack of 5 |
| 115-040518-00 | Mindray multifunction electrode (Paediatric) |
| 0651-30-77008 | Mindray multifunction electrode (Paediatric) Pack of 5 |
| Therapy Solutions | |
| 0651-30-77114 | External paddles kit, Adu/Paed, 1 conductive gel(250/tube) |
| 040-000545-00 | Therapy cable with test load (50 Ω) |
| Internal Paddles (Spoons) | |
| 115-018366-00 | Internal Defibrillator Paddle, Reusable, 1 inch |
| 115-018367-00 | Internal Defibrillator Paddle, Reusable, 2 inch |
| 115-018368-00 | Internal Defibrillator Paddle, Reusable, 3 inch |
| Carry cases & Mounting solutions | |
| 115-018610-00 | BeneHeart D3 Carry Case with Protection Shield |
| 115-066638-00 | Vehicle Mounting Kit |
| 115-013412-00 | Table mounting |
| 115-015825-00 | Trolley kit |
| Replacement Battery | |
| 115-007858-00 | Battery, Li-ion, (14.8V, 3000mAh) BeneHeart D3 |
| 115-049328-00 | Battery, Li-ion, (15.1V, 5600mAh) BeneHeart D3 Platinum only |
| Paper | |
| A30-000001- --- | 1 x Roll Thermal Recorder Paper |