

TE7 Max Ultrasound System

Transducer Family



SC5-1Ns (Convex)

Applications: Abdomen, Nerve, OB/GYN, Pediatrics, Urology, Vascular, Thoracic/Pleural
FOV: 100°
Convex radius: 60 mm
Biopsy kit: NGB-022, multi-angle, reusable



L11-3VNs (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NA
eSpatialNavi™: Enabled



C4-1s (Convex)

Applications: Abdomen, Cardiac, OB/GYN
FOV: 96°
Convex radius: 30 mm
Biopsy kit: NGB-036, multi-angle, reusable



L9-3s (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 43 mm
Biopsy kit: NGB-034, multi-angle, reusable



C11-3s (Micro-Convex)

Applications: Abdominal, Cardiac, Cranial, Nerve, Vascular
FOV: 100°
Convex radius: 15 mm
Biopsy kit: NGB-018, multi-angle, reusable



L20-5s (Linear)

Applications: MSK, Nerve, Small Parts, Superficial, Vascular
FOV: 28.5 mm
Biopsy kit: NA



L12-3VNs (Linear)

Applications: MSK, Nerve, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NA
eSpatialNavi™: Enabled



L14-5sp (Linear)

Applications: Intraoperative, MSK, Nerve, Superficial, Vascular
FOV: 25.3 mm
Biopsy kit: NA



L12-3RCs (Linear)

Applications: MSK, Nerve, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-043, multi-angle, reusable



L14-6Ns (Linear)

Applications: Abdominal, MSK, Cranial, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-007, multi-angle, reusable

**L14-6s (Linear)**

Applications: Pediatrics, MSK, Cranial, Nerve, Small Parts, Superficial, Vascular
FOV: 25 mm
Biopsy kit: NGB-016, multi-angle, reusable

**V11-3Ws (Endocavity)**

Applications: OB/GYN, Urology
FOV: 173°
Convex radius: 11 mm
Biopsy kit: NGB-004, single-angle, reusable

**SP5-1Ns (Phased)**

Applications: Abdominal, Cardiac, Cranial, Obstetrics, Thoracic/Pleural
FOV: 90°
Biopsy kit: NGB-011, multi-angle, reusable

**P7-3Ts (TEE)**

Applications: Transesophageal Echo
FOV: 90°
Biopsy kit: NA

**SP5-1s (Phased)**

Applications: Abdominal, Cardiac, Cranial, Thoracic/Pleural
FOV: 90°
Biopsy kit: NGB-011, multi-angle, reusable

**P8-3Ts (TEE)**

Applications: Transesophageal Echo
FOV: 90°
Biopsy kit: NA

**P7-3s (Phased)**

Applications: Abdominal, Cardiac, Cranial, Pediatric
FOV: 90°
Biopsy kit: NA

**P8-2Ts (TEE)**

Applications: Transesophageal Echo
FOV: 90°
Biopsy kit: NA

**P10-4s (Phased)**

Applications: Abdominal, Cardiac, Cranial
FOV: 90°
Biopsy kit: NA

**7LT4s (Linear)**

Applications: Abdominal, Intraoperative, MSK, Cranial, Pediatrics, Small Parts, Superficial, Vascular
FOV: 43 mm
Biopsy kit: NGB-010, multi-angle, reusable

Mindray North America

800 MacArthur Boulevard, Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 www.mindray.com

©2022 Mindray DS USA, Inc. Subject to change. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. 9/22 P/N: 0002-08-40539 Rev C

mindray
NORTH AMERICA