

# TE7 Ultrasound System

## Transducer Family



### C5-1s (Convex)

**Applications:** Adult Abdominal, Nerve, OB/GYN, Pediatrics, Urology, Vascular  
**FOV:** 101°  
**Convex radius:** 60 mm  
**Biopsy kit:** NGB-022, multi-angle, reusable



### L12-3RCs (Linear)

**Applications:** MSK, Nerve, Pediatrics, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NGB-043, multi-angle, reusable



### C5-2s (Convex)

**Applications:** Adult Abdominal, Nerve, OB/GYN, Pediatrics, Urology, Vascular  
**FOV:** 76°  
**Convex radius:** 51 mm  
**Biopsy kit:** NGB-015, multi-angle, reusable



### L12-4s (Linear)

**Applications:** MSK, Nerve, Pediatrics, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NGB-007, multi-angle, reusable



### C4-1s (Convex)

**Applications:** Abdomen, Cardiology, OB/GYN  
**FOV:** 96°  
**Convex radius:** 30 mm  
**Biopsy kit:** NGB-036, multi-angle, reusable; CIVCO 698-019, disposable



### L11-3VNs (Linear)

**Applications:** MSK, Nerve, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NA  
eSpatialNavi™: enabled



### C11-3s (Micro-Convex)

**Applications:** Neonatal Abdominal, Neonatal Cardiac, Neonatal Cranial, Nerve, Neurology, Pediatric Abdominal, Pediatric Cardiac, Vascular  
**FOV:** 100°  
**Convex radius:** 15 mm  
**Biopsy kit:** NGB-018, multi-angle, reusable



### L9-3s (Linear)

**Applications:** MSK, Nerve, Small Parts, Vascular  
**FOV:** 43 mm  
**Biopsy kit:** NGB-034, multi-angle, reusable



### L14-6Ns (Linear)

**Applications:** MSK, Nerve, Pediatrics, Small Parts, Superficial, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NGB-007, multi-angle, reusable



### 7L4s (Linear)

**Applications:** Adult Abdominal, MSK, Neonatal Cranial, Pediatrics, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NGB-007, multi-angle, reusable

**L7-3s (Linear)**

**Applications:** Adult Abdominal, MSK, Nerve, Pediatrics, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NGB-007, multi-angle, reusable

**L20-5s (Linear)**

**Applications:** MSK, Nerve, Small Parts, Superficial, Vascular  
**FOV:** 28.5 mm  
**Biopsy kit:** NA

**L14-6s (Linear)**

**Applications:** Pediatrics, MSK, Neonatal Cranial, Nerve, Small Parts, Superficial, Vascular  
**FOV:** 25 mm  
**Biopsy kit:** NGB-016, multi-angle, reusable

**L14-5sp (Linear)**

**Applications:** Intraoperative, MSK, Nerve, Superficial, Vascular  
**FOV:** 25.3 mm  
**Biopsy kit:** NA

**P10-4s (Sector)**

**Applications:** Neonatal Abdominal, Neonatal Cardiac, Neonatal Cranial  
**FOV:** 90°  
**Biopsy kit:** NA

**P7-3s (Sector)**

**Applications:** Abdominal, Cardiac, Cephalic, Pediatric  
**FOV:** 90°  
**Biopsy kit:** NA

**SP5-1s (Sector)**

**Applications:** Adult Abdominal, Adult Cardiac, Adult Cranial, Pediatric Cardiac  
**FOV:** 90°  
**Biopsy kit:** NGB-011, multi-angle, reusable

**P4-2s (Sector)**

**Applications:** Adult Abdominal, Adult Cardiac, Adult Cranial, Pediatric Cardiac, Pediatrics, Neonatal Cranial  
**FOV:** 90°  
**Biopsy kit:** NGB-011, multi-angle, reusable

**V11-3ws (Endocavity)**

**Applications:** OB/GYN, Urology  
**FOV:** 173°  
**Convex radius:** 11 mm  
**Biopsy kit:** NGB-004, single-angle, reusable

**P7-3Ts (TEE)**

**Applications:** Adult Transesophageal Echo  
**FOV:** 90°  
**Biopsy kit:** NA

**7LT4s (Linear)**

**Applications:** Adult Abdominal, Intraoperative, MSK, Neonatal Cranial, Pediatrics, Small Parts, Superficial, Vascular  
**FOV:** 43 mm  
**Biopsy kit:** NGB-010, multi-angle, reusable

**Mindray North America Innovation Center**

2100 Gold Street  
 San Jose, CA 95002

Tel: 800.288.2121 Support: 877.913.9663 [www.mindray.com](http://www.mindray.com)

©2018 Mindray DS USA, Inc. Subject to change. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. 10/18 P/N: 0002-08-40352 Rev B

**mindray**  
 healthcare within reach