

# TE X Ultrasound System

## Transducer Family



### SC6-1s (Convex)

**Applications:** Abdominal, Gynecology, Obstetrics, Vascular, Nerve, Musculoskeletal, Urology, Thoracic/Pleural, Small Parts, Pediatric  
**FOV:** 71°  
**Convex radius:** 60 mm  
**Biopsy guide:** NGB-022, multi-angle, reusable



### L12-3VNs (Linear)

**Applications:** Small Parts, Nerve, Musculoskeletal, Vascular, Thoracic/Pleural, Ophthalmic, Abdominal, Pediatric  
**FOV:** 38.2 mm  
**Biopsy guide:** NA



### C11-3s (Micro-Convex)

**Applications:** Abdominal, Vascular, Small Parts, Musculoskeletal, Pediatric, Cardiac, Cranial, Thoracic/Pleural  
**FOV:** 113°  
**Convex radius:** 15 mm  
**Biopsy guide:** NGB-018, multi-angle, reusable



### L9-3s (Linear)

**Applications:** Small Parts, Musculoskeletal, Vascular, Nerve, Thoracic/Pleural Abdominal, Pediatric, Obstetrics  
**FOV:** 43.9 mm  
**Biopsy guide:** NGB-034, multi-angle, reusable



### V11-3Hs (Endocavity)

**Applications:** Gynecology, Obstetrics, Urology  
**FOV:** 210°  
**Convex radius:** 11 mm  
**Biopsy guide:** NGB-025, single-angle, reusable



### L14-6Ns (Linear)

**Applications:** Small Parts, Nerve, Musculoskeletal, Vascular, Thoracic/Pleural, Ophthalmic, Abdominal, Pediatric  
**FOV:** 38.2 mm  
**Biopsy guide:** NGB-007, multi-angle, reusable



### L12-3RCs (Linear)

**Applications:** Small Parts, Nerve, Musculoskeletal, Vascular, Thoracic/Pleural, Ophthalmic, Abdominal, Pediatric  
**FOV:** 38.2 mm  
**Biopsy guide:** NGB-043, multi-angle, reusable; NGB-044, multi-depth, reusable



### L20-5s (Linear)

**Applications:** Musculoskeletal, Vascular, Nerve, Abdominal, Small Parts, Ophthalmic  
**FOV:** 28.7 mm  
**Biopsy guide:** NA



### **P8-2s (Phased)**

**Applications:** Pediatric Cardiac, Pediatric Abdominal, TCI, Neonatal Cardiac, Neonatal Abdomen, Neonatal Head  
**FOV:** 90°  
**Biopsy guide:** NA



### **P7-3Ts (TEE)**

**Applications:** Transesophageal Echo  
**FOV:** 90°  
**Biopsy guide:** NA



### **SP5-1s (Phased)**

**Applications:** Cardiac, Cephalic, Abdominal, Vascular, Obstetrics, Thoracic/Pleural  
**FOV:** 90°  
**Biopsy guide:** NGB-011, multi-angle, reusable



### **L16-4Hs (Intra-operative)**

**Applications:** Musculoskeletal, Nerve, Small Parts, Vascular, Pediatric, Intra-operative, Cephalic  
**FOV:** 25.4 mm  
**Biopsy guide:** NA



### **i3P Wireless (Phased)**

**Applications:** Cardiac, Abdominal, Lung, Cranial  
**FOV:** 90°  
**Biopsy guide:** NA

#### **Mindray North America**

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2024 Mindray DS USA, Inc. Subject to change. 01/24 P/N: 0002-08-40630 Rev 2.0

# **mindray**

healthcare within reach