



# Panorama ViewStation or WorkStation Default Settings

Hospital: \_\_\_\_\_

Department: \_\_\_\_\_

Clinical / Sales / Service Representative: \_\_\_\_\_

Date: \_\_\_\_\_

Software Version: \_\_\_\_\_

Required Signature: \_\_\_\_\_

\* WorkStation Option

		SYSTEM SETUP		
<b>Display</b>		Single	Dual	
<b>Number of Tiles per Display</b>		8	12	
<b>Laser Printer #1</b>		HP type _____		
<b>Laser Printer #2</b>		HP type _____		
PARAMETER COLORS		OPTIONS AVAILABLE		ACCOUNT SELECTION
	HR	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	Resp	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	T1	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	T2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	Delta Temp	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	SpO2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	SvO2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP1	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP3	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP4	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP5	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP6	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP7	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	IBP8	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	NIBP	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	T3	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	CO2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	O2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		
	N2O	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue		

PARAMETER COLORS		OPTIONS AVAILABLE	ACCOUNT SELECTION
	Agent	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	CO/CI	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	CCO/CCI	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	SVR	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	Tb	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	PAWP	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	SVRI	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	ScVO2	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	STAT CCO/CCI	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	PVC/min	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	Pause/min	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	
	ST	Green, Lt Green, Purple, Dk Red, Salmon, Orange, White, Yellow, Lt Blue, Lt Orange, Grey, Dk Blue	

PRINT SETUP		OPTIONS AVAILABLE		ACCOUNT SELECTION
	Laser Jet 1	Reports	Strip	
	Laser Jet 2	Reports	Strip	
VOLUME		OPTIONS AVAILABLE		ACCOUNT SELECTION
	System Alarms	Low	Med High	
	Physio Alarms			
	Orange	Low	Med High	
	Green	Low	Med High	
	Blue	Low	Med High	
SYSTEM ALARMS		OPTIONS AVAILABLE		ACCOUNT SELECTION
	Password Protection*	ON	OFF	
	Mute All	10 sec, 15 sec, 30 sec, 45 sec, 60 sec, 120 sec, Permanent		
	Tech. Sound	ON	OFF	
PASSWORDS		OPTIONS AVAILABLE		ACCOUNT SELECTION
	Alarm*	Maximum of 15 Characters		
	System	Maximum of 15 Characters		
DATE AND TIME		OPTIONS AVAILABLE		ACCOUNT SELECTION
	Date	YYYY-MM-DD; DD-MM-YYYY; MM-DD-YYYY:		
	Time	HH:MM(24H), HH:MM:SS(24H); HH:MM(12H); HH:MM:SS(12H),		
UNIT CHOICES		OPTIONS AVAILABLE		ACCOUNT SELECTION
	Trends Access Order	List, Graphic		
Demographics Order		Last Name, First Name, Bed, ID, Demo Off		
	Position 1	Last Name, First Name, Bed, ID, Demo Off		
	Position 2	Last Name, First Name, Bed, ID, Demo Off		
	Position 3	Last Name, First Name, Bed, ID, Demo Off		
	Position 4	Last Name, First Name, Bed, ID, Demo Off		
	Disable ECG	On	Off	
CLAN COMMUNICATION FROM PANORAMA CENTRAL STATION		Full	View Off	