

Datascope Accutorr® V



Refer to the Accutorr V Monitor Operation Manual for detailed user instructions.

Patient Info



Create a Patient ID
Change Patient Size

Patient Size



Select Patient Size

Set Alarms



Set Alarms

Start NIBP



Start NIBP
Measurement

Display



Toggle between Trend
and Waveform Display

Deflate



Stop NIBP
Measurement and
Deflate Cuff

Menu



System Setup

Interval



Set Interval

Silence



Silence Alarm

Print



Print
Stop Print

REV ID: 2001101170906							Adult	
Date	Time	Sys	Dia	Map	SpO2	PR	Temp	
01-01	23:06	132	80	86	98	73	98	
01-01	23:05	114	76	88	97	76	98	
01-01	23:04	121	75	101	98	79	98	
01-01	23:03	132	74	100	97	77	97	
01-01	23:02	125	67	87	98	79	98	

01-01-2011 17:11:38

Display for Trend Data, Pulse
Waveform, Messages and
Configuration Settings



Scan Trends, Select
and Execute Settings



Power

Mindray DS USA, Inc.
800 MacArthur Blvd., Mahwah, NJ 07430
Tel: 1.800.288.2121 or 201.995.8000
www.mindray.com

mindray[™]
NORTH AMERICA